
Sistema Telecomandato Ribaltabile ed Elevabile  
*Elevating Remote Controlled Tilting System*

# Helios

## Tavolo telecomandato elevabile con ribaltamento $\pm 90^\circ$ Brand new remote controlled tilting table $\pm 90^\circ$

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*It has been designed with a privilege to the aspects of lightness, compactness and reduced dimensions. Each table movement has been optimized to facilitate the use by the operator and to ensure maximum patient safety, preventing his movements during the examination.*

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*The focal distance (SID) is adjustable from 105 to 180 cm. Chest examinations are thus performed directly on the table with no need of an external bucky.*

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
*Unique in the world, Helios SFD is set to be upgraded into Helios DRF version. The modular mechanics allow for the replacement of the SFD/II with the flat panel holder group only by a few, rapid adjustments. A genial and innovative idea will allow everyone to purchase a direct to digital system now or in a near future.*

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*The X-ray tube motorized rotation of  $\pm 180^\circ$  allows the directing of the beam towards any external receptor devices (such as a stretcher). It is activated from remote or from the two control boards placed on SFD/DRF front or tube front.*

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*The tables are supplied with one or more retractable grids having different focal distance. It is possible to carry out examinations with the following modes:*



- Oscillating grid (SFD only)
- Fixed grid (SFD and DRF)
- Without grid, for pediatric or extremities (SFD and DRF).

*In the DRF version, the grids can be comfortably parked into two rails placed below the detector holder*

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*The tube motorized rotation is suitable to the execution of tomography and oblique projections without the mechanical bar support. The tomography (up to 5 tilting angles) is bidirectional and can be done in each position of the table and SFD/II group.*

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*Elevated patient coverage  
The large dimensions of the tabletop (250 x 70 cm) linked to the longitudinal excursion of the X-ray tube column and spot film device (or flat panel) group allows a patient to a total scan of 202 cm (SFD version) and 209 (DRF version). In this case, there is no need to reposition the patient on the table*

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*The table control console includes a touch screen panel on which all commands are activated and 4 joystick s which duplicate the table main controls (tilting, lateral displacement, column displacement, etc). Additionally all commands are duplicated on the control boards placed on SFD /DRF front and tube front.*

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#### Table with analogue images acquisition system

*CCD camera 512x512 or 1Kx1K of excellent quality which is connected to a high definition image intensifier from 9" to 12". In its basic version, the system comes with Last Image Hold function (LIH), recursive filter and one LCD 17" medical monitor. Optionally the system can be equipped with additional or larger monitors.*

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#### Table with digital images acquisition and processing device

*It is realized by a high resolution CCD camera 1024x1024 and by an image intensifier from 12" to 16". It is possible to carry out single spots or digital movies both in general diagnostics, gastrointestinal, orthopedic and peripheral angiography (DSA software package is optional). The included workstation is easy to use and offers a variety of acquisition, display and processing settings. The system is fully compliant to DICOM 3.0 protocol.*

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\* Optional

*Helios DRF supports a vast range of fluoroscopic and angiographic applications:*

- Gastrointestinal examinations
- Urogenital tract examinations
- Endoscopy
- Arthrography
- Venography
- Lymphography
- Myelography
- Step Angio
- Digital angiography \*
- Digital subtracted angiography \*

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\* Optional

*Helios DRF supports an extensive range of radiographic examinations:  
Skeletal radiography:*

- Skull
- Thorax
- Abdomen
- Vertebral column
- Pelvis
- Upper and lower extremities
- Tomography
- Long image display or stitching\*

\* Optional



**Nuovo flat panel dinamico 43x43 cm**  
*new dynamic flat panel 43x43 cm*

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*The dynamic flat panel with active area 43x43cm is suitable for the examination of any anatomic region and provides almost 50% more coverage than a 16" image intensifier and about 20% than a 35x43 CR cassette.*

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*The Flat Panel technology delivers sharp and accurate images to support the diagnosis and is not affected by geometric distortion which is typical of image intensifiers.*

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Fluoroscopia in tempo reale fino a 30 imm/secondo  
*Real time fluoroscopy up to 30 fps*

Radiografia e fluoroscopia a campo pieno  
*Full field radiography and fluoroscopy*

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
## "2 in 1" efficiency

With Helios DRF the time needed to perform mixed fluoroscopy and radiography examinations is considerably reduced, using now the same system with direct to digital technique.

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## Reduce time

The use of cassettes and film developers tends to disuse. With Helios DRF the images are accessed in real time, shortening examination times for the operator and the patients.

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## Optimized connectivity

Thanks to its DICOM 3.0 capabilities, the system can be seamlessly integrated with existing HIS/RIS and PACS. If the objective is a filmless environment, it is possible to create CD's and DVD's directly from console with the integrated CD/DVD writer, or alternatively to send to print only the most significant images.

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The stitching function consists of an automatic reconstruction of an X-ray image starting from a series of images acquired at fixed frequency during the scanning of the patient. The image is reconstructed, keeping all original pixels, and can be viewed on the monitor, processed, printed or sent to the network. This function is usually used for the scan of spine and lower arts. Using the "stitching" command, the system automatically generates a sequence of acquisitions according to the number of frames selected by the operator. (2 for 60 cm, 3 for 90 cm, 4 for 120 cm length). Afterwards, the software links the images to get a single one shown in approximately 5 seconds.

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The minimum height from ground of the patient tabletop is 43 cm. The access is made easier for any kind of patient.

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#### DRF + Wi-Fi


Helios DRF can be combined with ASTER DR 1 FLAT Wi-Fi system which includes a ceiling tube stand, a vertical bucky holder and a 35x43 cm cassette-sized wireless flat panel. With the use of a wireless portable detector and a ceiling tube stand the system achieves maximum flexibility making imaging possible virtually everywhere in the room, as conducting examinations on patients in wheelchairs or stretchers. Adding a vertical bucky stand would complete the spectrum of upright examinations such as lateral projections. A single acquisition console is able to capture images coming from both dynamic flat panel and from the wireless one.

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#### SFD + Wi-Fi

*Even Helios SFD table is ready to be integrated with ASTER DR 1 FLAT Wi-Fi system. The 35x43cm flat panel is inserted in the spot film device like a conventional cassette and the images are captured and sent Wi-Fi to a dedicated console. When choosing for stitching mode (optional feature), images reconstructions will be done exactly like DRF version. HELIOS SFD is definitely the only system capable to make the stitching through a sport film device.*

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
*The Tomosynthesis is a technique which allows the reconstruction of volumetric images starting from a finite number of bi-dimensional projections taken at different tube angles. This is given by the separation of the acquisition process from the visualization one. Among the advantages with respect to CT scan technique, it is well accepted by those patients affected by severe claustrophobia. Tomosynthesis has several areas of application including:*

- Weight bearing views
- Knee pain
- Spontaneous Osteonecrosis
- Bone on bone
- Worn cartilage

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*BETA RF Generators ensure an elevated quality standard for the acquisition of radiological and fluoroscopic information coming from TV chain (Helios SFD Analog) or from digital systems (Helios SFD Digital and Helios DRF). Main characteristics are:*

- Power supply from 40 to 100 kW, frequency up to 400 KHz
- 1024 APR can be programmed in fluoroscopy + 768 in radiography (standard membrane console).
- AEC interface board for up to 4 ion chambers \*
- Supply for an additional x-ray tube \*
- Anode high speed starter \*
- Touch screen console with 20.000 APR programmable by the operator \*

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