

Intelligent Blocker Ice 1000



Ice 1000

Melody of Intelligence

– Enhanced intelligence provides the maximum comfort in the complex process of blocking lenses

The *Ice 1000* is a highest-class automatic blocker ever.

processing a lens. It works in perfect combination with the NIDEK lens edger families, the Lex 1000, and the Me 1200.





Measurement accuracy is at its highest level with the newly designed optical systems of the latest technology. Lenses with weak power can be measured properly and progressive lenses can be detected automatically. Errors caused by prismatic effect are eliminated.







Single Vision Mode

Multi-focal Mode

Progressive Mode

Lens Clamping Mechanism

The *Ice1000* features a unique lens clamp which secures the lens automatically. It enables the blocking of super-hydrophobic lenses without risk of axis shift.



Open position for Integrated Shape Imager



Closed position for lens blocking (multiposition lens stage is motorized)



Set the lens on the stage



Clamp keeps lens in a steady position while blocking

Integrated Shape Imager with hole edit function

With lens measurement table, The Ice 1000 accurately and automatically reads the outline shape of a demo lens.

Hole positions are set on the touch panel using a demo lens. Setting hole positions is easy, as images are enlarged on the display.

Step 1 Measure the lens shape.


Select the hole tab and set the demo lens on the table. The *Ice 1000* automatically and accurately digitizes the lens shape without the need of painting the edge by a marker.


(Some particular demo lenses may require painting on the edge for automatic shape reading.)



Step 2 Measure the hole position automatically.

The *Ice1000* can automatically read the circumference, hole position and horizontal axis of the demo lens from the image taken with a COMS camera. Processing data can be made extremely easily.






Step 3 Check the data and the job is complete.

Input data will be displayed for review. Data can be saved as a pattern.

■ High-performance editing

The *Ice 1000* has an array of editing functions that can support the *Me 1200* complex operations such as hole creation, Advanced Shape Editor, Design Cut, Partial Grooving and Facet. It is well suited for managing the edited data for processing on the *Me 1200*. Thumbnail images are shown in the data management display which is helpful when searching for missing data.








Layout

Blocking

Shape Editor








Design Cut

Partial Grooving

Facet


■ Information Bar

Even first-time users can easily operate the *Ice1000* using the Information Bar which provides useful "next-step" information.



■ Color-coded Lens Identification

To aid the lens edging process, the display of the *Ice 1000* utilizes a color-coding system adapted from navigational light colors, to identify right and left lenses: green represents the right lens and red represents the left lens. Red and green cups are included. You can avoid processing the wrong lens by following the color-prompted indicator on the display.



This unit has two RS-232C ports and a LAN port. This unit can drive multiple devices

and operate as a server for small to medium sized labs.



Rx Data Management

The *Ice1000* can store and serve Rx data which contains holes, Design Cut, Partial Grooving and Facet. Data viewer is convenient for finding stored shape data.



USB Flash Drive Port

Rx data is stored to USB flash drive, which is unlimited based on the size of the USB flash drive device and or PC hard drive the *Lex 1000* to the *Ice 1000* is also possible. Data management is easy and flexible.



Rotating Cup Adapter

The cup adapter is turned right-side up, making accurate cup insertion easy.

■ Large Touch Panel Display

The *Ice 1000* features a large color display making it extremely user-friendly. Even first time users can easily operate the unit, following the "next step" instructions of the Information Bar.

Color Camera

The unit is equipped with a color camera with high accuracy and auto-brightness control. It works well with all tinted lenses, gradient tints and polarized lenses.

■ Built-in Barcode Scanner (optional)

The Ice 1000's integrated barcode scanner saves lab space and streamlines the lens finishing process.



You can choose a side panel color to best suit your practice. A beautiful Lavender color (standard) or stylish Silver Mica* color are available.

*option only available in some countries





Lavender (standard)


Silver Mica (option)



System feasibility

■ Connection to lens edgers

The *Ice 1000* can be connected to the *Lex 1000*, *Me 1200*, and *SE-9090* series. The *Ice 1000* is compatible with the VCA protocols.



Ice 1000 / Ice 1000NT Specifications

Model	Ice 1000	Ice 1000NT
Tracer	Built-in	None
Lens size	Dia. 80 mm or less	←
Layout span	FPD: 30.0 to 99.5 mm	
	PD (or 1/2 PD): 30.0 to 99.5 mm (15.0 to 49.75 mm)	←
	Height of the optical center: 0 to ±15.0 mm	-
	Size adjustment: 0 to ±10.0 mm	
Item to be entered	FPD	
	PD (or 1/2 PD)	
	Cylinder axis	
	EP (eye point height of progressive lens)	
	Lens size	
	Lens material (Plastic, Plastic with high refractive index,	←
	Glass, Polycarbonate, Acrylic, TRIVEX)	
	Frame type (Metal, Celluloid, Nylor, Rimless)	
	Grinding mode selection	
	CYL (+/- switching)	
	Job code	
Lens measuring mode	Single vision mode: Full auto / Mark detection	
	Multi focal mode: Segment detection	
	Progressive mode: Print mark / Print mark angle / Point mark	←
	Manual mode	
	Demo lens block mode	
Shape imager function	Measurement range: 65.0 x 50.0 mm (±1.5 mm)	
	Hole position: 0.01 mm increments	←
	Hole diameter: Ø0.5 to 10.0 mm (0.01 mm increments)	
Blocking method	Auto	←
Power supply	AC 100 to 120, AC 200 to 240 V, 50 / 60 Hz	←
Power consumption	100 VA	←
Dimensions / Mass	325 (W) x 507 (D) x 345 (H) mm / 24.0 kg	325 (W) x 507 (D) x 345 (H) mm / 18.0 kg
	12.8 (W) x 20.0 (D) x 13.6 (H) " / 52.9 lbs.	12.8 (W) x 20.0 (D) x 13.6 (H) " / 39.7 lbs.
Standard accessories	Power cord, Spare fuse, Interface cable, Stylus pen,	Power cord, Spare fuse, Interface cable, Stylus pen,
	Pattern setting unit, Standard frame, Standard pattern,	Accessory case, Blower, Frame change holder,
	Accessory case, Blower, Frame change holder, USB flash drive,	USB flash drive, Lens measurement table
	Lens measurement table	
Optional accessories	Barcode scanner	←

Specifications and design are subject to change without notice.



HEAD OFFICE

34-14 Maehama, Hiroishi Gamagori, Aichi 443-0038, Japan Telephone:+81-533-67-6611 Facsimile:+81-533-67-6610 URL: http://www.nidek.co.jp

[Manufacturer]

TOKYO OFFICE (International Div.)

3F Sumitomo Fudosan Hongo Bldg., 3-22-5 Hongo, Bunkyo-ku, Tokyo 113-0033, Japan

Telephone: +81-3-5844-2641 Facsimile :+81-3-5844-2642 URL: http://www.nidek.com NIDEK INC.

47651 Westinghouse Drive Fremont, CA 94539, U.S.A. Telephone:+1-510-226-5700

:+1-800-223-9044 (US only) Facsimile : +1-510-226-5750 URL : http://usa.nidek.com

NIDEK S.A.

Europarc 13, rue Auguste Perret 94042 Créteil, France

Telephone: +33-1-49 80 97 97 Facsimile :+33-1-49 80 32 08 URL: http://www.nidek.fr **NIDEK TECHNOLOGIES Srl**

Via dell'Artigianato, 6 / A 35020 Albignasego (Padova), Italy Telephone: +39 049 8629200 / 8626399 Facsimile : +39 049 8626824

URL: http://www.nidektechnologies.it

