ERES  
MEDICAL®

# Spirobank II Smart

Available on the  
App Store

Bluetooth®  
SMART

**MIR**  
MEDICAL INTERNATIONAL RESEARCH

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
**Recommendations**

**Key recommendations**

- FVC Phase 1: Breathing at rest (Tidal Volume)
- FVC Phase 2: Rapid Maximum Inspiration
- FVC Phase 3: Respiratory Pause
- FVC Phase 4: Start of Forced Expiration
- FVC Phase 5: Peak expiratory flow
- FVC Phase 6: Maximum Expiration
- FVC Phase 7: Maximum Forced expiration
- Spirometry Quality

**Virtual Assistance (V.A.) control**

- FVC Phase 2: Inflow/hold rapid maximum inspiration
- FVC Phase 3: Excessive respiratory pause
- FVC Phase 4: Start of forced expiration (no leak detected)

**Age limit to carry out spirometry**

Spirometry can be performed from age 5 years (in some cases also 3-4 years). Manual reference and standard recommendations for performing spirometry please refer to the international ATS/ERS guidelines.

**Correct posture to carry out test**

It is preferable to perform it in a standing position, using a chair with armrests but without wheels, for two reasons. Firstly because most of reference values (normal or abnormal) were established in a standing position, secondly for safety reasons, because forced expiration reduces venous return, with the possibility of syncope or of dizziness in even more easily fainting. If the patient complains of dizziness, it is recommended to stop the spirometry and wait for a few minutes before repeating the spirometry. If the problem persists it is preferable to reschedule the session.

**Suggestions on how to carry out a vital capacity test (V.C)**

If the patient complains of dizziness, it is recommended stop the forced spirometry test and perform a vital capacity (VC) measurement waiting to repeat the forced spirometry test if possible. If no problem persists it is preferable to reschedule the session.

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