
INSTRUCTIONS FOR USE

AtmosAir

Mattress Replacement System

AtmosAir Fit
AtmosAit Plus

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WARNING

To avoid injury, always read this *Instructions for Use* and accompanied documents before using the product.

Mandatory to read the *Instructions for Use*.

® and ™ are trademarks belonging to the Arjo group of companies

© Arjo 2019.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

ARJO HEREBY DISCLAIMS ALL EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE ARJO PRODUCT AS DESCRIBED IN THIS PUBLICATION. UNDER NO CIRCUMSTANCES SHALL ARJO BE LIABLE FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES AND EXPENSES, INCLUDING DAMAGES OR INJURY TO PERSON OR PROPERTY, DUE IN WHOLE OR IN PART TO THE USE OF THE PRODUCT OTHER THAN THOSE FOR WHICH DISCLAIMER OF WARRANTY OR LIMITATION OF LIABILITY IS EXPRESSLY PROHIBITED BY SPECIFIC, APPLICABLE LAW. NO PERSON HAS THE AUTHORITY TO BIND ARJO TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH IN THIS PARAGRAPH.

Descriptions or specifications in Arjo printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Information in this publication may be subject to change at any time. Contact Arjo for updates.

IMPORTANT INFORMATION FOR USERS

In order for Arjo products to perform properly, Arjo recommends the following conditions. Failure to comply with these conditions will void any applicable warranties.

- Use this product only in accordance with these instructions and applicable product labeling.
- Assembly, operations, adjustments, extensions, modifications, technical maintenance or repairs must be performed only by qualified personnel authorized by Arjo. Contact Arjo for information regarding maintenance and repair.

Specific indications, contraindications, warnings, precautions and safety information exist for Arjo's therapeutic support systems. It is important for users to read and familiarize themselves with these instructions and to consult the treating physician prior to patient placement and product use. Individual patient conditions may vary.

Intentionally left blank

Table Of Contents

- Introduction.....2**
 - Indications 2
 - Contraindications..... 2
 - Risks and Precautions..... 2
 - Safety Information 3
 - Serious Incident 4
- Preparation for Use5**
- Mattress Installation.....6**
 - Installing Side Bolsters 6
 - Patient Placement and Nursing Care 7
 - CPR..... 7
 - Incontinence / Drainage 7
 - General Operation..... 7
- Care and Cleaning8**
 - Fire Barrier 8
 - Cover Cleaning Options 8
- Preventive Maintenance Schedule.....9**
 - Daily Cleaning 9
 - Inspection / System Check-Out..... 9
- Troubleshooting10**
- Parts Diagram – AtmosAir Fit MRS.....11**
- Parts Diagram – AtmosAir Plus MRS12**
- Replacement Parts13**
- Specifications14**
- Symbols Used.....15**
- Customer Contact Information.....15**

Introduction

CAUTION: *It is important that you carefully read and review these instructions with your health care professional. If you have any questions of a medical nature, please contact your health care professional.*

It is recommended that all sections of this user manual be read prior to product use. Carefully review the **Indications, Contraindications, Risks and Precautions** and **Safety Information** prior to placing a patient on any *AtmosAir™ Fit Mattress Replacement System (MRS)* and *AtmosAir™ Plus Mattress Replacement System (MRS)*.

Caregivers should review this information with the patient and the patient's family and/or legal guardian. Save this User Guide in an easily accessible location for quick reference.

The *AtmosAir Fit MRS* and *AtmosAir Plus MRS* are non-powered pressure redistribution mattress replacement systems which uses *Self Adjusting Technology (SAT)* to provide pressure redistribution therapy. The systems are designed for patients weighing up to 1000 lb (454 kg), with expandability features from 36 in (91 cm) without bolsters to 42 in to 48 in (107 cm to 122 cm) wide with bolsters on the *AtmosAir Fit MRS* and 34 in (86 cm) without bolsters to 41 in to 48 in (104 cm to 122 cm) wide with bolsters on the *AtmosAir Plus MRS*.

Indications

The *AtmosAir Fit MRS* and *AtmosAir Plus MRS* are indicated for the prevention and treatment of skin breakdown.

Contraindications

- unstable vertebral fracture
- cervical and skeletal traction

Risks and Precautions

Transfer – Standard precautions should be taken during patient transfer.

Side Rails and Restraints – **WARNING:** Use or non-use of restraints, including side rails, can be critical to patient safety. Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of side rails or other restraints. **See related Safety Information.**

Patient Migration – Specialty surfaces have different shear and support characteristics than conventional surfaces and may increase the risk of patient movement, sinking and / or migration into hazardous positions of entrapment and / or inadvertent bed exit.

Monitor patients frequently to guard against patient entrapment.

Safety Information

Skin Care – Monitor skin conditions regularly and consider adjunct or alternative therapies for high acuity patients. Give extra attention to any possible pressure points and locations where moisture or incontinence may occur or collect. Early intervention may be essential to preventing skin breakdown.

Patient Weight – The maximum patient weight for these devices is 1000 lb (454 kg). In addition, consult the specifications for the bed frame being used. Additional weight limitations may apply.

Patient Entrance / Exit – Caregiver should always aid patient in exiting the bed. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency.

Brakes – Caster brakes should always be locked once the bed is in position. Verify wheels are locked before any patient transfer to or from the bed.

Bed Height – To minimize risk of falls or injury, the bed should always be in the lowest practical position when the patient is unattended.

Bed Frame – Always use a standard bariatric healthcare bed frame with these mattresses, with any safeguards or protocols that may be appropriate. Bed frame and side rails (if used) must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body. In the U.S. it is recommended that bed and side rails (if used) comply with the FDA's [Hospital Bed System Dimensional and Assessment Guidance To Reduce Entrapment](#).

Head of Bed Elevation – Keep head of bed as low as possible to help prevent patient migration.

Side Rails / Patient Restraints – Whether and how to use side rails or restraints is a decision that should be based on each patient's needs and should be made by the patient and the patient's family, physician and caregivers, with facility protocols in mind. Caregivers should assess risks and benefits of side rail / restraint use (including entrapment and patient falls from bed) in conjunction with individual patient needs and should discuss use or non-use with patient and / or family. This includes assessment of the bed occupant and the combination of bed frame, side rail and mattress (or mattresses where overlays are used). Risk assessment should be repeated if the bed frame, mattress, side rail or condition of patient changes. Consider not only the clinical and other needs of the patient but also the risks of fatal or serious injury from falling out of bed and from patient entrapment in or around the side rails, restraints or other accessories. In the US, for a description of entrapment hazards, a description of at-risk patients and guidance to further entrapment risks, refer to the FDA's [Hospital Bed System Dimensional and Assessment Guidance To Reduce Entrapment](#). Outside the US, consult the local Competent Authority or Government Agency for Medical Device Safety for specific local guidance. Consult a caregiver and carefully consider the use of bolsters, positioning aids or floor pads, especially with confused, restless or agitated patients. It is recommended that side rails (if used) be locked in the full upright position when the patient is unattended. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency. **Monitor patients frequently to guard against patient entrapment.**

CAUTION: *When selecting a bariatric mattress, ensure the distance between top of side rails (if used) and top of mattress (without compression) is at least 8.66" (220 mm) to help prevent inadvertent bed exit or falls. Consider individual patient size, position (relative to the top of the side rail) and patient condition in assessing fall risk.*

CPR – Level the bed. Lower side rails and initiate CPR per facility protocols. Consider use of backboard if indicated. After CPR remove backboard, if used, raise siderails and reconfigure bed and accessories as in initial placement.

No Smoking in Bed – Smoking in bed can be dangerous. To avoid the risk of fire, smoking in bed should never be allowed.

General Protocols – Follow all applicable safety rules and institution protocols concerning patient and caregiver safety.

Fire Barrier Contact – Care should be taken with regards to the fire barrier when the cover is removed. It is recommended that Personal Protective Equipment (PPE) is used during handling of the mattress with the fire barrier exposed to avoid contact with loose particles. To ensure optimal product performance, avoid over-handling of the fire barrier and handle with care. If the fire barrier becomes damaged or soiled it should be replaced.

End of Life Disposal –

- Fabric material used on the mattresses or any other textiles, polymers or plastic materials etc. should be sorted as combustible waste.
- Mattresses at the end of life should be disposed of as waste according to the national or local requirements which may be landfill or combustion.
- Pump units have electrical and electronic components should be disassembled and recycled per Waste of Electrical and Electronic Equipment (WEEE) or in accordance with local or national regulation.

Serious Incident

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

Preparation for Use

For information concerning the bed frame refer to the manufacturers' User Guide.

1. Open shipping container(s).

Do not use sharp instruments to open boxes. Damage to mattress could result.

2. Remove the *AtmosAir* Fit MRS or *AtmosAir* Plus MRS from plastic protective cover.

The mattress cover may appear wrinkled when unpacked. To remove wrinkles, allow mattress up to 24 hours to acclimate; see Troubleshooting table for more information. Wrinkles will not affect inflation or function, so MRS may be used immediately if needed.

3. Check mattress surface for tears or cracking; do not use if tears or cracks are present.
4. If re-installing MRS onto a new frame or for a new patient, check mattress surface for staining and soiling; clean and / or disinfect as required (see **Care and Cleaning**).
5. Level bed and lock brakes.
6. Remove existing mattress from bed frame.

Mattress Installation

1. Position mattress on bed frame with logo facing up and product information tags at foot end of bed.
2. Ensure mattress is properly positioned with no gaps between mattress and bed frame or side rails.

Always use a standard bariatric health care bed frame with safeguards or protocols that may be appropriate. Frame and side rails must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body.


Installing Side Bolsters

After expanding the bariatric bed, add the bolsters using the following procedure:

1. Place the bolster in the gap between the side rails and the mattress with the zippers to the outside and the marked foot end at the foot end of the bed.
2. Connect the four female buckles on the underside of the left and right side bolsters to the four male buckles on the underside of the mattress.

AtmosAir Plus MRS side bolsters are side-specific. The left bolster (identified by an inward-facing foot) must go on the patient left side of the bed. Same applies to the right side bolster.

3. Zip the mattress to the mattress inserts.

When bolsters are not in use, tuck bolster sleeves under mattress.

Patient Placement and Nursing Care

It is recommended that all sections of this manual be read prior to patient placement and nursing care. Carefully review the **Contraindications**, **Safety Information** and **Risks and Precautions** sections prior to placing a patient on any mattress.

1. Transfer patient following all applicable safety rules and institution protocols.
2. Center patient side-to-side and head-to-foot on mattress surface.
3. Ensure all sections of the mattress fully support the patient.

Mattress handles are for mattress transport only.

CPR

1. Level bed.
2. Lower or remove side rails as necessary.
3. Begin CPR as per facility protocol. Consider use of backboard if indicated.
4. After CPR is performed:
 - Remove backboard if used.
 - Raise or install side rail as necessary.
 - Reconfigure bed and accessories as in initial placement.

Skin Care

- Remove excess moisture and keep skin dry and clean.
- Check patient's skin regularly, particularly in areas where incontinence and drainage occur.
- Ensure linens under patient are not wrinkled.

Incontinence / Drainage

- Use moisture-impermeable underpads for incontinent patients.
- Wipe surface clean and replace bed linens as required (see **Care and Cleaning** if needed).

General Operation

Avoid contact of sharp instruments with the mattress. Punctures, cuts and tears may prevent proper inflation and air pressure maintenance.

Care and Cleaning

The following processes is recommended, but should be adapted to comply with local institution protocols. If you are uncertain, you should seek advice from your local Infection Control Specialist.

The *AtmosAir* MRS should be routinely decontaminated between patients and at regular intervals while in use.

Do not use Phenol-based solutions or abrasive compounds or pads on the cover during decontamination process as these will damage the surface coating. Do not boil or autoclave the cover.

Covers should be wiped down immediately after it has been exposed to fluids or spills.

Fire Barrier

Care should be taken with regards to the fire barrier when the cover is removed. It is recommended that Personal Protective Equipment (PPE) is used during handling of the mattress with the fire barrier exposed to avoid contact with loose particles. To ensure optimal product performance, avoid overhandling of the fire barrier and handle with care. Fire barrier should be replaced if soiled or damaged.

Cover Cleaning Options

Sewn *AtmosAir* MRS (Top cover can not be removed from base)

1. Remove or push bed linens to center of mattress.

The sewn *AtmosAir* MRS cover is cleaned by wipe down methods only. Do not launder as damage to the base may occur.

2. Wipe and rinse any soiling from the mattress surface and base. Use 1000 ppm chlorine or 70% alcohol.
3. Rinse with a wet towel to remove chemicals.
4. After wipe down, dry surface with towel.
5. Ensure bed linens are refitted and not wrinkled under patient.
6. Clean pump and tubes (if needed) by wiping with a damp cloth.

RF Welded *AtmosAir* MRS (Detachable top cover)

1. Push bed linens to center of mattress to wipe down, remove to launder top cover or wipedown as described below.

The RF Welded *AtmosAir* MRS has a detachable top cover that can be laundered. However the base is cleaned by wipe down methods only. Do not launder base as damage may occur.

2. Wipe and rinse any soiling from the mattress surface and base. Use 1000 ppm chlorine or 70% alcohol.
3. Rinse with a wet towel to remove chemicals.
4. After wipe down, dry surface with towel.
5. Ensure bed linens are refitted and not wrinkled under patient.

Laundrying Detachable Top Cover

1. Unzip the top cover from the base for laundrying.
Do not launder base as damage may occur.
2. Recommended wash temperature for top cover is 60° C (140° F) for 15 minutes.
3. Maximum wash temperature is 95° C (203° F) for 15 minutes.
4. Tumble dry at 60° C (140° F) or air dry.
5. Maximum drying temperature 80° C (176° F).

Preventive Maintenance Schedule

Preventive maintenance for the *AtmosAir* Fit MRS or *AtmosAir* Plus MRS consists of regular cleaning (see **Care and Cleaning**) and an overall system check-out to be performed at the intervals described below.

All components must be cleaned, disinfected and inspected after each patient's use and before use by a new patient. Always use standard precautions, treating all used equipment as contaminated. Institutions should follow local protocols for cleaning and disinfection.

Daily Cleaning

The cover should be wiped daily with a mild soap and water solution. After wipe down, dry surface with a towel.

Inspection / System Check-Out

Check each of the following before placing the mattress with a new patient:

1. Check mattress surface for tears or cracking; do not use if tears or cracks are present.
2. Ensure mattress is free of stains and is not overly faded.

Troubleshooting

It is recommended that all sections of this manual be reviewed before troubleshooting any *AtmosAir* Fit MRS or *AtmosAir* Plus MRS.

Do not attempt troubleshooting outside this manual or where the remedy recommends to contact an Arjo service representative. Any unauthorized service, modification, alteration or misuse may lead to serious injury and / or product damage and will void all applicable warranties.

| SYMPTOM | POSSIBLE CAUSE | SOLUTION |
|--|---|---|
| Mattress too firm upon arrival | Difference in altitude not sufficient to open valves. | Apply weight to mattress to open valves. |
| Mattress cover too wrinkled upon removal from shipping container | Internal components have not acclimated to environment. This does not effect inflation or function. | Let mattress acclimate for 24 hours. If problem continues, contact Arjo for assistance. |
| Mattress does not inflate or is not firm. | Tubing not connected properly. | Check tubing inside mattress for loose connectors. |
| | Tubing kinked. | Check tubing inside mattress for kinks. |
| | Tubing disconnected. | Check tubing inside mattress for possible disconnect. |
| | Holes in or damage to SAT system. | Check SAT system for holes or damage or contact Arjo for assistance. |




Parts Diagram – AtmosAir Fit MRS

All hoses were removed to improve diagram readability.

The fire barrier is a sleeve that fits over the foam.

The **SAT Cover**, **SAT Cells** and **SAT Interliner** are all part of the integrated **SAT System** and cannot be ordered separately. See **Replacement Parts** for a complete list of **SAT Systems**.


Parts Diagram – AtmosAir Plus MRS

All hoses were removed to improve diagram readability.




The fire barrier is a sleeve that fits over the foam.

The **SAT Cover**, **SAT Cells** and **SAT Interliner** are all part of the integrated **SAT System** and cannot be ordered separately. See **Replacement Parts** for a complete list of **SAT Systems**.

Sewn Mattress
Cover and Base
Cover not detachable

Replacement Parts

Do not attempt troubleshooting, maintenance or parts replacement outside this manual or where the remedy recommends contacting an Arjo service representative. Any unauthorized service, modification, alteration or misuse may lead to serious injury and / or product damage and will void all applicable warranties.

Replaceable *AtmosAir* Fit MRS components are listed below. For more information such as pricing or additional spare parts that are not on this list, please contact your local Arjo representative.

AtmosAir Fit MRS

| PART TYPE..... | CATALOG NUMBER |
|---|-----------------------|
| <i>AtmosAir</i> Fit MRS (includes 1 mattress and 2 side bolsters) | 312493 |

AtmosAir Fit MRS Replacement Parts

| PART TYPE..... | CATALOG NUMBER |
|-----------------------|-----------------------|
| Side Bolster | 312489 |
| Cover | AFCLRVG35080TXS |
| Mattress..... | 312492 |

AtmosAir Plus MRS

| PART TYPE..... | CATALOG NUMBER |
|--|-----------------------|
| <i>AtmosAir</i> Plus MRS (includes 1 mattress, 1 right side bolster, 1 left side bolster, and 1 foot bolster)..... | APMRLVG34079TFS |

AtmosAir Plus MRS Replacement Parts

| PART TYPE..... | CATALOG NUMBER |
|--------------------------------------|-----------------------|
| Base Mattress Replacement | APRRLVG34079TFS |
| Base Mattress Replacement Cover..... | APCRLVG34079TXS |
| Right Side Bolster..... | APBRLVG07079TFD |
| Right Side Bolster Cover | APCRLVG07079TXD |
| Left Side Bolster | APBRLVG07079TFL |
| Left Side Bolster Cover..... | APCRLVG07079TXL |
| Foot Bolster | APBRLVG48509TXF |
| Foot Bolster Cover..... | APCRLVG48509TXS |

Specifications*

Maximum Weight Capacity 1000 lb (454 kg)

AtmosAir Fit MRS Mattress:

Mattress Weight 56 lb (25 kg)

Mattress Weight with bolsters 80 lb (36 kg)

Mattress Length 80 in (203 cm)

Mattress Width with bolsters 42 in or 48 in (107 cm or 122 cm)

Mattress Width without bolsters 35 in (91 cm)

Mattress Height 7 in (18 cm)

AtmosAir Fit MRS Bolsters:

Bolsters 80 in x 6 in (203 cm x 15 cm)

Bolster Weight 12 lb (5 kg) each

Bolster Height 7 in (18 cm)

AtmosAir Plus MRS Mattress:

Mattress Weight 55 lb (25 kg)

Mattress Weight with bolsters 87 lb (39 kg)

Mattress Length 79.5 in (202 cm)

Mattress Width with bolsters 41 in or 48 in (104 cm or 122 cm)

Mattress Width without bolsters 34 in (86 cm)

AtmosAir Plus MRS Bolsters:

Bolsters (L x W x H) 79.5 in x 7 in x 9 in (202 cm x 18 cm x 23 cm)

Bolster Weight 16 lb (7 kg) each

***Specifications subject to change without notice.**

Symbols Used

Important Operational Information

Foot End

Warning of possible hazard to system, patient or staff

Consult User Guide

Manufacturer

Safe Working Load

Wash at 60°C for 15 minutes,
Max 95°C for 15 minutes

Serial Number

The operator must read this document (this instruction for use) before use. Note: the symbol is blue on the product label.

Do Not Iron

Use Solution diluted to 1000 ppm of available chlorine

Attention – See User Guide

CE marking indicating conformity with European Community harmonised legislation

Wipe down Only

Tumble Dry at 60°C (140°F),
Max 80°C (176°F)

Do not use Phenol-based cleaning solutions

Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745

Customer Contact Information

For questions regarding this product, supplies, maintenance or additional information about Arjo products and service, please contact Arjo or an Arjo authorized representative or visit: www.arjo.com.

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Intentionally left blank

Intentionally left blank

AUSTRALIA
Arjo Australia Pty Ltd
78, Forsyth Street
O'Connor
AU-6163 Western Australia
Tel: +61 89337 4111
Free: +1 800 072 040
Fax: +61 89337 9077

BELGIQUE / BELGIË
Arjo NV/SA
Evenbroekveld 16
BE-9420 ERPE-MERE
Tél/Tel: +32 (0) 53 60 73 80
Fax: +32 (0) 53 60 73 81
E-mail: info.belgium@arjo.com

BRASIL
Arjo Brasil Equipamentos Médicos
Ltda
Rua Marina Ciufuli Zanfelice, 329
PB02 Galpão - Lapa
São Paulo – SP – Brasil
CEP: 05040-000
Phone: 55-11-3588-5088
E-mail: vendas.latam@arjo.com
E-mail: servicios.latam@arjo.com

CANADA
Arjo Canada Inc.
90 Matheson Boulevard West
Suite 300
CA-MISSISSAUGA, ON, L5R 3R3
Tel/Tél: +1 905 238 7880
Free: +1 800 665 4831 Institutional
Free: +1 800 868 0441 Home Care
Fax: +1 905 238 7881
E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA
Arjo Czech Republic s.r.o.
Na Strži 1702/65
140 00 Praha
Czech Republic
Phone No: +420225092307
e-mail: info.cz@arjo.com

DANMARK
Arjo A/S
Vassingerødvej 52
DK-3540 LYNGBE
Tel: +45 49 13 84 86
Fax: +45 49 13 84 87
E-mail:
dk_kundeservice@arjo.com

DEUTSCHLAND
Arjo GmbH
Peter-Sander-Strasse 10
DE-55252 MAINZ-KASTEL
Tel: +49 (0) 6134 186 0
Fax: +49 (0) 6134 186 160
E-mail: info-de@arjo.com

ESPAÑA
Arjo Ibérica S.L.
Parque Empresarial Rivas Futura, C/
Marie Curie 5
Edificio Alfa Planta 6 oficina 6.1.-62
ES-28521 Rivas Vacia, MADRID
Tel: +34 93 583 11 20
Fax: +34 93 583 11 22
E-mail: info.es@arjo.com

FRANCE
Arjo SAS
2 Avenue Alcide de Gasperi
CS 70133
FR-59436 RONCQ CEDEX
Tél: +33 (0) 3 20 28 13 13
Fax: +33 (0) 3 20 28 13 14
E-mail: info.france@arjo.com

HONG KONG
Arjo Hong Kong Limited
Room 411-414, 4/F,
Manhattan Centre,
8 Kwai Cheong Road,
Kwai Chung, N.T.,
HONG KONG
Tel: +852 2960 7600
Fax: +852 2960 1711

ITALIA
Arjo Italia S.p.A.
Via Giacomo Peroni 400-402
IT-00131 ROMA
Tel: +39 (0) 6 87426211
Fax: +39 (0) 6 87426222
E-mail: Italy.promo@arjo.com

MIDDLE EAST
Arjo Middle East FZ-LLC
Office 908, 9th Floor,
HQ Building, North Tower,
Dubai Science Park,
Al Barsha South
P.O Box 11488, Dubai,
United Arab Emirates
Direct +971 487 48053
Fax +971 487 48072
Email: Info.ME@arjo.com

NEDERLAND
Arjo BV
Biezenwei 21
4004 MB TIEL
Postbus 6116
4000 HC TIEL
Tel: +31 (0) 344 64 08 00
Fax: +31 (0) 344 64 08 85
E-mail: info.nl@arjo.com

NEW ZEALAND
Arjo Ltd
34 Vestey Drive
Mount Wellington
NZ-AUCKLAND 1060
Tel: +64 (0) 9 573 5344
Free Call: 0800 000 151
Fax: +64 (0) 9 573 5384
E-mail: nz.info@Arjo.com

NORGE
Arjo Norway AS
Olaf Helsetts vei 5
N-0694 OSLO
Tel: +47 22 08 00 50
Faks: +47 22 08 00 51
E-mail: no.kundeservice@arjo.com

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Address page - REV 24: 04/2019

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjo.com

arjo

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Arjo Inc.