

Low Temperature Sterilizers

PLASMASTER 40 / 80 / 140

The innovative Belimed H₂O₂ low-temperature plasma sterilizers ensure safe and reliable sterilization of thermolabile materials and instruments. The advantages include short processing times, no outgassing time, no danger to people or the environment, gentle treatment of valuable instruments, dedicated and competent service – all from the same supplier.



The PLASMASTER sterilizer 140 is available as one-door (see picture) or two-door version (PLASMASTER 140D). The two-door version is completely integrable in the CSSD.

The procedure is quick and safe

Thanks to the highly effective, safe and environmentally friendly Belimed PLASMASTER, the sterilized material is available for use immediately after sterilization. Outgassing time is not necessary because no toxic substances are used. The H₂O₂ used is converted to natural water vapor and oxygen.

Easy operation and full data traceability

The user-friendly touch panel and clearly spoken status and warning notifications ensure easy operation. Thanks to the standard integrated USB port and printer, all cycle data is recorded and printed. Belimed PLASMASTER sterilizers can also be integrated with the optional Belimed ICS 8535 batch documentation system.

Easy installation, intelligent design

The built-in rollers allow the sterilizer to be easily positioned for use. Only a single phase connection with 230 volts at 50 hertz is required and the sterilizer is ready to operate. Optimized usable chamber volume supports the economical usage of the sterilizer. Via an Ethernet interface technical support and remote diagnostic can be provided.

Belimed
Infection Control

System solutions for cleaning, disinfection and sterilization in the healthcare, pharmaceutical and laboratory sectors

Belimed Austria: +43 3155 40699 0, Belgium: +32 10 42 02 40, China: +86 21 513 709 98, France: +33 4 37 41 63 03, Germany: +49 8631 9896 0, Hungary +36 1 318 86 97, Netherlands: +31 10 286 17 50, Slovenia: +386 1 7866 010, Switzerland: 0848 55 88 11, United Kingdom: +44 1403 738 811, USA: +1 843 216 7424, other countries: +41 41 449 78 88, www.belimed.com

Technical data

	PLASMASTER 40	PLASMASTER 80	PLASMASTER 140	PLASMASTER 140D
Chamber				
Chamber shape	square	circular	square	square
Chamber dimensions H x W x D (mm)	210 x 330 x 640	Ø 400 x 650	425 x 425 x 790	425 x 425 x 790
Chamber volume (l)	44	80	142	142
Useable chamber volume (l)	40	71	120	120
Installation				
External dimensions H x W x D (mm)	1600 x 760 x 825	1675 x 752 x 835	1671 x 831 x 1012	1685 x 1055 x 1075
Weight (kg)	approx. 350	approx. 430	approx. 600	approx. 600 (one-door) approx. 650 (two-door)
Doors	1, manual	1, manual	1, manual	2, automatic
Operation	5.6" full color touch panel with graphic text display and voice activation	5.6" full color touch panel with graphic text display and voice activation	8.4" full color touch panel with graphic text display and voice activation	10.4" full color touch panel with graphic text display and voice activation
Documentation	integrated batch documentation printer, optional connection to the Belimed Infection Control Software ICS 8535	integrated batch documentation printer, optional connection to the Belimed Infection Control Software ICS 8535	integrated batch documentation printer, optional connection to the Belimed Infection Control Software ICS 8535	integrated batch documentation printer, optional connection to the Belimed Infection Control Software ICS 8535
Consumption of H₂O₂ per batch (ml)	4.0 ml	7.0 ml	10.0 ml	10.0 ml

Standard programs

Procedure	Standard program	PLASMASTER 40	PLASMASTER 80	PLASMASTER 140	PLASMASTER 140D
Vacuum phase 1	Initial evacuation of the chamber (min)	4–6	6–7	8–20	8–10
Diffusion phase 1	Automatic injection and diffusion phase (min)	11–13	11–12	21–22	21–22
Vacuum phase 2 and plasma phase 1	Repeated evacuation of the chamber (min)	1–1,5	1–2	1–2	1–2
Diffusion phase 2	Automatic injection and diffusion phase (min)	11–13	11–12	21–22	21–22
Plasma phase 2	Repeated evacuation of the chamber (min)	1–1,5	1–2	2–3	2–3
Degassing and drying phase	Adjustment of the sterilization chamber to room pressure and repeated evacuation of the chamber (min)	4–5	4–5	7–8	7–8
	Total cycling time (min)	35–40	42–45	65–70	65–70



PLASMASTER 40 PLASMASTER 80 PLASMASTER 140 PLASMASTER 140D

Standards and guidelines met

1. European Medical Products Guideline 93/42/EEC
2. ISO 9001:2000 Quality Management Systems – Design, Development and Manufacture
3. ISO 13485:2003 Medical Products – Quality Management Systems – Design, Development and Manufacture CE 0120
4. Can be validated in accordance with ISO 14937

Belimed
Infection Control

System solutions for cleaning, disinfection and sterilization in the healthcare, pharmaceutical and laboratory sectors

Belimed Austria: +43 3155 40699 0, Belgium: +32 10 42 02 40, China: +86 21 513 709 98, France: +33 4 37 41 63 03, Germany: +49 8631 9896 0, Hungary +36 1 318 86 97, Netherlands: +31 10 286 17 50, Slovenia: +386 1 7866 010, Switzerland: 0848 55 88 11, United Kingdom: +44 1403 738 811, USA: +1 843 216 7424, other countries: +41 41 449 78 88, www.belimed.com