Dati tecnici / Technical data

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

MINIMUM HARDWARE AND SOFTWARE REQUIREMENTS

PC Desktop:

Processor Intel Pentium Dual Core 2.00 GHz - 2 GB RAM - Firewire 1394A board OHCI 1.1 compatible - Video board 512 MB RAM (not shared) resolution 1280x1024 pixels

Software:

Operating system Microsoft Windows XP Home SPK3, Windows XP Professional SPK3, Windows Vista 32 bit Home premium, Windows 7 Home and Pro – 32 bit and 64 bit.

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
THE NEXT GENERATION
OF NON-MYDRIATIC FUNDUS CAMERA

Costruzione Strumenti Oftalmici
Via degli Stagnacci 12/E, 50018 Scandicci
Firenze, ITALY
Tel. ++39 055.722191, Fax. ++39 055.721557
www.csoitalia.it

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HIGH RESOLUTION ACQUISITION
COBRA can acquire a very detailed picture with the integrated digital color-camera.

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COBRA

Cobra: the next Generation of non-mydriatric fundus camera

COBRA is an innovative non-midriatic digital fundus camera that combines all the needed functions for a fast retinal screening. By using a new optical system, COBRA provides high quality retinal pictures. Designed with a very hergonomic shape, COBRA provides a clear and detailed vision of the entire ocular fundus, with a real 60° horizontal field of view. Pictures are acquired with a minimum flash exposition, thus by minimizing the examination time for the patient. The high resolution CCD sensor (5 MP in the HD version) is both used for the patient alignment (with an IR source) and for the retinal picture acquisition (with a white LED flash and IR). COBRA is connected to the PC by IEEE1394 protocol (Firewire type A or B), providing a very high speed data transfer. Patients data are saved in a database, using the Phoenix software. Stand-alone or Network station are both possible. DICOM feature is integrated in Phoenix software.

Phoenix Software imaging and archiving features

image processing advanced functions
Drawing features
Measurement tools
Mosaic function
Meibomiam Glands Dinsfunction (MGD) analysis tool
AVr evaluation for heart attack risk
Advanced printout functions (comparison visible-IR and OD-OS at the same time)
Wavelengths splitter summary printout on the same page.
Integration with ERG test performed by Retimax

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

RETIMAX OVERLAP TOOL

COBRA provides the new RETIMAX plus Fundus tool, that adds new attractive feature of the Multifocal ERG test, by combining Multifocal ERG and photographic image of the retinal Fundus. This new feature helps in the detection and following the progression of a macular or other limited retinal area, with specific and accurate indication of the functionality of each retinal ar

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
INFRARED PICTURE ACQUISITION
InfraRED picture is acquired by using an LED InfraRED light source.

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
MOSAIC FUNCTION
By using the manual acquisition, Cobra can acquire pictures from different retinal zones with the possibility to show a wider retinal vision thanks to the MOSAIC function.

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WAVELENGTHS SPLITTER
The wavelengths splitter shows in a screen: in the upper part the original picture, the infrared one and the red-free; in the bottom part the choroidal component, the vascular and the optic fiber.

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AVr EVALUATION TOOL
The new Avr tool provides an advanced and complete screening of the cardio-vascular health of the patient, by performing an analysis of the health of arteriorials and veins in the ocular fundus. An index called Avr (Arteriolar-To-Venular ratio) provides a risk percentage in the period of 3 years.

MGD ANALYSIS TOOL
COBRA include san advanced tool for Meibomiam Glands Disfunction (MGD) analysis. This tool, by using an easy color picture acquisition, and without any additional accessory, provides an evaluation of the glands loss criticity (5 steps scale).

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