clear valley separator of the platelets and RBC peaks. Then there was no clear separation of platelets and RBC based on size. This is often due to platelet aggregation in cells and that cells often have large platelets that merge with their small RBC in volume.


Figure 1: Photo shows normal erythrocytes, a monocyte (M1 side) and lymphocyte.

Moreover the erythrocytes were normochromic and normochromic, the platelets were of normal size. The blood smear had erythrocytes of normal size without any increase in immature cells. Fewer erythrocytes are smaller than normal and show much less central pallor but much less central pallor than that usually seen in the normal dog.

There was a very clear valley on the blood RBC histogram showing great separation of the lymphocytes and granulocytes. This shows that the valley pattern suggested the advanced lymphocyte and granulocyte counts should be accurate. The relative lymphocyte count increased and the relative granulocyte count decreased. The relative granulocyte count increased and the relative lymphocyte count decreased. The absolute number was within reference values. The manual monocyte count and automated red cell cell did not match well.



Figure 2: Note 2 normal neutrophils, erythrocytes and a lymphocyte (lower left).

There were only 16 monocytes seen on the smear at differential leukocyte count as the smear at differential leukocyte count as the smear did not reflect this. The absolute monocyte count for manual and automated monocyte counts to vary. The absolute monocyte count is not the same as monocyte count. Monocytes are often not evenly distributed on a blood smear. It is common that a manual monocyte count per 100 RBCs is an automated monocyte count per 1000 RBCs. The manual monocyte count and automated red cell cell did not match well.

Figure 3: 3 pale neutrophils have red shaped granules.

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