

PRIME VERSIONS WITH AIR DRYER








50/25 PRIME M








Code: PR-0F300-040-M

Description: Oil free compressors, FAD at 5 bar 225 l/min, tank capacity 40 l. with internal surface treatment against corrosion, equipped with the new membrane air dryer able to warranty till -20°C dew point value

This compressor is available in 10 BAR VERSION also! visit the "10 bar oil-free compressors section" for more details!

Model	 Dental chairs	 Air dryer	 Air Capacity	 FAD @ 5bar	 Motor	 Volt. Freq.	 Receiver
			m ³ /h l/min cfm	m ³ /h l/min cfm	HP KW	V Hz	l GAL
50/25 PRIME M	4	yes	- 360 12,6	- 225 7,9	3 2,2	230-400 50	40 10,57

Model	 Noise Level	 Duty cycle	 Max pressure	 Size	 Gross weight
	dB(A)		bar PSI	mm	Kg lb
50/25 PRIME M	66	70%	7 101	710X410X810	58 127

