


## Disposable Microtome Blades

- › S35
- › R35
- › N35
- › N35HR
- › A35
- › A22
- › S22
- › C35
- › S35L
- › S35LL
- › High Profile



pfm medical ag offers a broad range of FEATHER® disposable microtome blades for paraffin and cryostat sectioning for the whole range of individual applications in routine practice, research and industry.

[www.pfmmedical.com](http://www.pfmmedical.com)



pfm medical ag offers a broad range of **FEATHER® disposable microtome blades** for paraffin and cryostat sectioning for the whole range of individual applications in routine practice, research and industry.

### S35



#### Use

The S35 disposable microtome blade is a standard blade for routine practice and biopsies.

#### Details

- ▶ With excellent sectioning quality from the beginning
- ▶ Very long service life
- ▶ Short trimming time

#### Technical Data


- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

#### Ordering Information

**Type S35**

**REF 207500000**

### R35



#### Use

The R35 disposable microtome blade is extra suitable for routine ribbon sectioning and for hard tissues.

#### Details

- ▶ Specially hardened cutting edge
- ▶ Short trimming time

#### Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

#### Ordering Information

**Type R35**

**REF 207500005**

## N35



### Use

The N35 disposable microtome blade is particularly suitable for hard tissues.

### Details

- ▶ Longer service life
- ▶ Suitable for all routine sectioning, also for soft tissue and biopsies
- ▶ Also suitable for ribbon sectioning after trimming
- ▶ Low compression

### Technical Data


- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

### Ordering Information

#### Type N35

REF 207500006

## N35HR



### Use

The N35HR disposable microtome blade is particularly suitable for sectioning hard and difficult tissues.

### Details

- ▶ Suitable for routine and ribbon sectioning
- ▶ Low compression
- ▶ Longer service life
- ▶ Short trimming time
- ▶ New coating
- ▶ Specially finished cutting edge

### Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

### Ordering Information

#### Type N35HR

REF 207500014

**A35**



**Use**

The A35 disposable microtome blade allows the production of very thin sections and is suitable for hard and soft tissues.

**Details**

- ▶ Very long service life (also for hard tissue)
- ▶ Short trimming time with corresponding clearance angle
- ▶ Very sharp

**Technical Data**

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

**Ordering Information**

**Type A35**

**REF 207500011**

**A22**



**Use**

The A22 disposable microtome blade allows the production of very thin sections in routine use and research and is suitable for fibrous tissues (biopsies/fabrics).

**Details**

- ▶ Long service life
- ▶ Short trimming time with corresponding clearance angle
- ▶ Very sharp
- ▶ Also suitable for ribbon sectioning after trimming

**Technical Data**

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 22°
- ▶ **Dispenser:** 50 pcs.

**Ordering Information**

**Type A22**

**REF 207500002**

**S22**



**Use**

The S22 disposable microtome blade allows the production of extra-thin sections and is suitable for soft tissues/fabrics.

**Details**

- ▶ Long service life
- ▶ Short trimming time with corresponding clearance angle

**Technical Data**

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 22°
- ▶ **Dispenser:** 50 pcs.

**Ordering Information**

**Type S22**

**REF 207500001**

**C35**



**Use**

The C35 disposable microtome blade is specifically designed for cryostat sectioning.

**Details**

- ▶ Carbon steel blade
- ▶ Short trimming time

**Technical Data**

- ▶ **Material:** carbon steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 20 pcs.

**Ordering Information**

**Type C35**

**REF 207500003**

## S35L



### Use

The S35L disposable microtome blade with a length of 120 mm is particularly suitable for large tissue blocks and routine sectioning.

### Details

- ▶ Short trimming time

### Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 120 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 20 pcs.

### Ordering Information

**Type S35L**

**REF 207500004**

## S35LL



### Use

The extra-long S35LL disposable microtome blade with a length of 180 mm is suitable for large tissue blocks and routine sectioning.

### Details

- ▶ Short trimming time

### Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 180 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 20 pcs.

### Ordering Information

**Type S35LL**

**REF 207500009**

## High Profile



### Use

The High Profile disposable microtome blade is suitable for all paraffin embedded tissue blocks in routine applications.

### Details

- ▶ Short trimming time

### Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 14 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

### Ordering Information

**Type High Profile**

REF 207500010

## Blade Depot



### Use

The blade depot provides suitable storage for FEATHER® disposable microtome blades.

### Details

- ▶ Storage capacity: approx. 200 blades
- ▶ Accessorie: protective cover blade deposit

### Manufacturer

pfm medical ag  
Wankelstraße 60  
50996 Köln, Germany

### Ordering Information

**Blade Depot**












REF 030000

**Protective cover blade deposit**

REF 030100

Ordering Information

Disposable Microtome Blades

Type	Colour coding	REF	Material	Length	Height	Blade angle	PU
S35	 dark blue	207500000	stainless steel	80 mm	8 mm	35°	50 pcs.
R35	 black	207500005	stainless steel	80 mm	8 mm	35°	50 pcs.
N35	 light green	207500006	stainless steel	80 mm	8 mm	35°	50 pcs.
N35HR	 dark green	207500014	stainless steel	80 mm	8 mm	35°	50 pcs.
A35	 light blue	207500011	stainless steel	80 mm	8 mm	35°	50 pcs.
A22	 pink	207500002	stainless steel	80 mm	8 mm	22°	50 pcs.
S22	 grey	207500001	stainless steel	80 mm	8 mm	22°	50 pcs.
C35	 white	207500003	carbon steel	80 mm	8 mm	35°	20 pcs.
S35L	 blue	207500004	stainless steel	120 mm	8 mm	35°	20 pcs.
S35LL	 green	207500009	stainless steel	180 mm	8 mm	35°	20 pcs.
High Profile	 brown	207500010	stainless steel	80 mm	14 mm	35°	50 pcs.

Blade Depot

REF	Capacity	PU
030000	200 blades	1 pc.

Protective cover blade deposit

REF	PU
030100	1 pc.

Contact

Should you have any questions our Customer Solutions Team is happy to advise you.

service@pfmmedical.com  
 T +49 (0)2236 9641-220  
 F +49 (0)2236 9641-51

www.pfmmedical.com - Your information source for the products of pfm medical ag.



pfm medical ag  
 Wankelstraße 60  
 50996 Köln, Germany

Exclusive Distributor  
 for Europe

Certified by  
 DIN EN ISO 13485

Manufacturer  
 FEATHER® Safety Razor Co. Ltd.



Registered Trademark in Japan  
 Products are subject to alterations