


Ventilator Support System

For Safe and Versatile Fixation of Ventilator Circuits



- Sterilizable for maximum patient protection against cross-contamination
- Height adjustable vertical post for varied positioning requirements
- Ergonomic joint handle for rapid set-up
- Versatile twin tube holder accepts 3 different tubing diameters
- Distal arm allows for tube orientation to prevent stress on patients mask
- Universal Clamp attaches to a vertical or horizontal post
- OR Table Rail Clamp available

DIMENSIONS

Vertical Post: 18.3" (46.5cm) x 5/8" (1.6cm)
 Mid Segment Arm: 17.3" (44cm)
 Distal Segment Arm: 8" (20cm)
 Tube Holder Capacity: 22mm, 16mm, 10mm

ORDERING INFORMATION

REF. NO.	DESCRIPTION
300-VSA01	Ventilator Support Arm with Universal Clamp (300-UC)
300-VSA02	Ventilator Support Arm with OR Table Rail Clamp (300-RC)
301-VSA	Ventilator Support Arm (no Rail Clamp)
300-UC	Universal Clamp
300-RC	OR Table Rail Clamp

Mediflex® Surgical Products
 250 Gibbs Road
 Islandia, NY 11749 USA

631.582.6424
info@mediflex.com
www.mediflex.com

MOUNTING CLAMPS

Universal Clamp

Ref. No.
300-UC



Mounts to vertical or horizontal posts
 (bed rails, IV poles, etc.)

OR Table Rail Clamp

Ref. No. 300-RC



Mounts to surgical table
 rail or rail segments