

Bookler[®] Endoscopic FlexArm[™] Plus

HOLDING & POSITIONING DEVICES



- Versatile endoscope and instrument holding device
- Expanded utility of Bookler[®] and Bookwalter[™] systems*
- Multiple configurations to accommodate all surgical practices
- Individual components offered to maximize flexibility

Complete system (Ref. No. 99190) in use with Quick-Grip Tip (10mm) (Ref. No. 72162). Order Separately.
Scope not included.

STRONGER • SMOOTHER • EASIER TO USE


Bookler® Endoscopic FlexArm™ Plus

The Bookler FlexArm Plus System provides extreme versatility as a surgical holding and positioning device. An assortment of Holding Tips are offered to accommodate any laparoscope or laparoscopic instrument. Individual system components are available to maximize the flexibility of the Bookler FlexArm Plus System.

COMPLETE BOOKLER® ENDOSCOPIC SYSTEMS

REF. NO.	DESCRIPTION
99190	Bookler FlexArm Plus with Quick-Connect Tip-12"(30cm) Post
99105	Bookler FlexArm Plus with fixed Stainless Steel Tip-12"(30cm) Post
99106	Bookler FlexArm Plus with fixed Delrin Tip-12"(30cm) Post
99191	Bookler FlexArm Plus with Quick-Connect Tip-18"(45cm) Post
99095	Bookler FlexArm Plus with fixed Stainless Steel Tip-18"(45cm) Post
99096	Bookler FlexArm Plus with fixed Delrin Tip-18"(45cm) Post

FlexArm Devices are available with a Quick-Connect (A) Hexagonal Fitting used with our interchangeable tips or can be ordered with permanently attached Stainless Steel Tip (B) or Delrin Tip (C).







Note: All Systems include a Post Coupler (Ref. No. 72013) and a 12" (30cm) or 18" (45 cm) Field Post

INDIVIDUAL COMPONENTS

BOOKLER FLEXARM PLUS (ARM ONLY)



◀ 99103	Bookler FlexArm Plus with Quick-Connect Hex Fitting (A)
99051	Bookler FlexArm Plus with fixed Stainless Steel Tip (B)
99029	Bookler FlexArm Plus with fixed Delrin Tip (C)

FIELD POSTS	POST COUPLING	VERTICAL EXTENSION KIT
 REF. NO. 72011 Standard 12" (30cm) REF. NO. 72009 Long 18" (45cm)	 REF. NO. 72013 Post Coupler 1" x 1" (2.5 x 2.5cm)	 REF. NO. 72193 12" Vertical Extension Bar with Post Coupler (72013)
	 REF. NO. 72013-C* Quick Lock Post Coupler 1" x 1" (2.5 x 2.5cm) <small>*not for sale in the U.S.</small>	



HOLDING TIPS

ADJUSTABLE HOLDING TIPS

 REF. NO. 69706 Stainless Steel Tip Holding Capacity: 3.5mm - 15mm	 REF. NO. 69707 Delrin Tip Holding Capacity: 3.0mm - 15mm	 REF. NO. 69706-L XL Stainless Steel Tip Holding Capacity: 5mm - 31.75mm
---	--	---

AUTOMATIC HOLDING TIPS

Safely and quickly secures 5mm & 10mm endoscopes. Atraumatic jaws will not damage scopes.

REF. NO. 72162 Quick-Grip Tip (10mm)		REF. NO. 72163 Quick-Grip Tip (5mm)	
--	---	---	---

*Bookler is a registered trademark of Flexbar Machine Corporation. Bookwaller is a registered trademark of Symmetry Surgical Inc.