

## HOLDING TIPS AND ACCESSORIES

FlexArm Devices are available with a Quick-Connect (A) Hexagonal Fitting used with our interchangeable tips or can be ordered with permanently attached Stainless Steel Tip (B) or Delrin Tip (C).



A - Quick-Connect Hexagonal Fitting



B - Permanently attached Stainless Steel Tip



C - Permanently attached Delrin Tip

### ADJUSTABLE HOLDING TIPS



REF. NO. 69706  
Stainless Steel Tip



REF. NO. 69707  
Delrin Tip





REF. NO. 69706-L  
XL Stainless Steel Tip

### AUTOMATIC HOLDING TIPS

REF. NO. 72162  
Quick-Grip Tip (10mm)



REF. NO. 72163  
Quick-Grip Tip (5mm)



Mediflex Holders and Positioners with a Quick-Connect distal end allows any Mediflex Tip equipped with a Hex Fitting to attach to any arm.

### STERILIZATION TRAYS



REF. NO. 69760  
Multi-Use Lidded Tray,  
Perforated, Aluminum  
Dimensions: 22.75" x 10.75" x 4.75"  
(57.7cm x 26.9cm x 11.9cm)



REF. NO. 69761  
Multi-Use Lidded Tray,  
Perforated, Aluminum  
Dimensions: 26.75" x 11" x 6"  
(67.8cm x 27.9cm x 15.2cm)



REF. NO. 69762  
Multi-Use Lidded Tray,  
Perforated, Aluminum  
Dimensions: 31" x 11" x 6"  
(78.7cm x 27.9cm x 15.2cm)

# Mediflex®

SURGICAL PRODUCTS

Mediflex Surgical Products,  
a Division of Flexbar Machine Corp.  
250 Gibbs Road  
Islandia, NY 11749 USA


Phone: 631.582.6424  
Toll- Free: 800.764.0474  
Fax: 631.582.4520  
Email: sales@mediflex.com

[mediflex.com](http://mediflex.com)

Made in U.S.A.

# Mediflex®

SURGICAL PRODUCTS



## FlexArm™ Plus

HOLDING & POSITIONING DEVICES

CE

Visit [mediflex.com](http://mediflex.com) for product specifications

FAP-0116

STRONGER • SMOOTHER • EASIER TO USE



## FlexArm™ Plus Holding & Positioning Devices

For hands-free, maximum strength instrument and scope holding which permits intra-operative adjustment, unlimited positioning and versatile instrument holding capability.

Mediflex, innovator of the original FlexArm Surgical Holding and Positioning Arms is pleased to introduce its latest generation technology. The new FlexArm Plus provides several engineering enhancements which result in a stronger, smoother and easier to use device.

### ENHANCED STRENGTH

New design results in a stronger arm with the ability to make adjustments when arm is tensioned.

### ERGONOMIC HANDLE DESIGN

Completely new handle design permits comfortable, quicker and easier tensioning.

### NO NEED TO RESET


By eliminating the reset procedure, it's now easier to clean and process which saves time as well as ensures optimal performance, every time.

### EASY SET-UP


New design permits an easier, quicker set-up. Fewer handle turns are required to go from loose to rigid.

### INCREASED WEAR RESISTANCE

Multiple engineering enhancements make the new arm design more robust with smoother actuation which in turn, extends the life of the product.



## MONOLITHIC FLEXARM™ PLUS



The Monolithic FlexArm Plus features a one-piece design with a built in OR Rail Clamp which is easy to set-up and simple to use. The ergonomic drop-down handle automatically aligns parallel to the arm to prevent interference.

REF. NO.	DESCRIPTION
99045	Monolithic FlexArm Plus with fixed Stainless Steel Tip
99045-D	Monolithic FlexArm Plus with fixed Delrin Tip
99045-QC	Monolithic FlexArm Plus with QC Hex Fitting
99045-QG10	Monolithic FlexArm Plus with 10mm Quick-Grip Tip
99045-QG5	Monolithic FlexArm Plus with 5mm Quick-Grip Tip
99702	Long Reach Monolithic FlexArm Plus with QC Hex Fitting*
99701	Nathanson Long Reach Monolithic FlexArm System with 5mm Hooks (see chart below)
99720	Nathanson Long Reach Monolithic FlexArm System with 6.5mm Hooks (see chart below)

\* Long Reach model approx. 4" (10.1cm) longer horizontal reach than 99045-QC system.

## NATHANSON FLEXARM™ PLUS SYSTEMS

For Bariatric and GERD procedures, the FlexArm Plus coupled with a Nathanson Liver Retractor permits static, atraumatic retraction of the liver without retractor sagging or a human retractor-holder fatiguing.

REF. NO.	DESCRIPTION	Ref. No. 99701 5mm System contains (1 ea)	Ref. No. 99720 6.5mm System contains (1 ea)
99702	Long Reach Monolithic FlexArm Plus	•	•
69703	Small 5mm Nathanson Hook	•	
69704	Med. 5mm Nathanson Hook	•	
69705	Large 5mm Nathanson Hook	•	
69735	Small 6.5mm Nathanson Hook		•
69736	Med. 6.5mm Nathanson Hook		•
69737	Large 6.5mm Nathanson Hook		•
69706	Stainless Steel Tip	•	•
69707	Delrin Tip	•	•

Note: Additional Nathanson Hook configurations available



## FLEXARM™ PLUS SINGLE ARMS & SYSTEMS

An assortment of FlexArm Plus System configurations are available for hands-free instrument and/or scope holding and stability during any surgical procedure which can reduce staff requirements and procedure time while increasing OR efficiencies.



REF. NO.	DESCRIPTION
99054-QC	Single Arm System with Quick-Connect Hex Fitting
99054	Single Arm System with fixed Stainless Steel Tip
99054-D	Single Arm System with fixed Delrin Tip
99055-QC	Twin Arm System with Quick-Connect Hex Fittings
99055	Twin Arm System with fixed Stainless Steel Tips
99055-D	Twin Arm System with fixed Delrin Tips
99055-DS	Twin Arm System with one fixed Stainless Steel Tip & one fixed Delrin Tip
99054-QCLR	Single Arm System, Long Reach Post with Quick-Connect Hex Fitting
99054-LR	Single Arm System, Long Reach with fixed Stainless Steel Tip
99054-DLR	Single Arm System, Long Reach with fixed Delrin Tip
99055-QCLR	Twin Arm, Long Reach Post w/Quick-Connect Hex Fittings
99055-LR	Twin Arm System, Long Reach, with fixed Stainless Steel Tips
99055-DLR	Twin Arm System, Long Reach, with fixed Delrin Tips
99055-DSLRL	Twin Arm System, Long Reach, with one fixed Stainless Steel Tip & one fixed Delrin Tip
99056-QC	Single Arm only with Quick-Connect Hex Fitting
99056	Single Arm only with fixed Stainless Steel Tip
99056-D	Single Arm only with fixed Delrin Tip
99056-QC-BS01	Single Arm only with Distal Swivel Option & Quick-Connect Hex Fitting

Note: All systems include an OR Rail Clamp.

## FLEXARM™ PLUS COMBO SYSTEM

The FlexArm Plus – StrongArm System features the absolute rigidity of a StrongArm device with the flexible rigidity of a FlexArm holder and positioner. The combination is ideal for procedures which require organ retraction while also holding a scope.

REF. NO.	DESCRIPTION
99040	StrongArm FlexArm Plus Combo System with Quick-Connect Hex Fittings and a heavy duty OR Rail Clamp



Visit [mediflex.com](http://mediflex.com) for product specifications