


# Thyroid Radiofrequency Ablation

Minimally Invasive Treatment for Thyroid Nodules



# STARmed Electrodes

- The sharpest electrodes that need no incision
- Clear visualisation of electrode tip on ultrasound
- Easy & smooth electrode control for "moving-shot" technique
- Complete range of electrode types & specifications



# Generator



- Advanced generator system
- Dedicated algorithm to prevent charring
- Impedance sound guide for better control of procedure


# About **STAR**med

## Global Leader in Thyroid RFA

- Largest number of published clinical studies performed with STARmed devices
- Partnership with Key Opinion Leaders for Thyroid RFA procedure training
- New device and clinical procedure development based on customer feedback
- Continuous professional educational activities held around the world

## Globally Recognised Treatment Option

- Pioneered Thyroid RFA into international guidelines (Austria, Italy, Korea, U.K., U.S.A.)
- Improves patient quality of life\*




# 65+

# CLINICAL ARTICLES

More than 65 Thyroid RFA articles have been published with STARmed electrodes. These articles prove the success of this treatment by the reduction in thyroid nodule size, symptom and cosmetic scores, and quality of life.

## \* Ordering Number Electrodes

Model Name / Specification	star Fixed	VIVA II	star Bipolar	star Uni_Faced	star Injectable
Gauge	18/19 G	18 G	17 G	17 G	18/19 G
Length	7-10 cm	7-10 cm	7-10 cm	7-10 cm	7-10 cm
Active Tip	3.8, 5, 7, 10, 15 mm	5-30 mm	12, 14, 15, 16 mm	5, 7, 10, 15 mm	7, 10, 15 mm



### Manufacturer



STARmed Co., Ltd.  
#B-14th Floor, Daebang Triplaon  
158 Haneulmaeul-ro, Ilsandong-gu  
Goyang-si Gyeonggi-do  
10355 Korea

### Distributor



inomed Medizintechnik GmbH  
Im Hausgruen 29  
79312 Emmendingen (GERMANY)  
Tel. +49 7641 9414-0 · Fax +49 7641 9414-94  
info@inomed.com · www.inomed.com