

ASKIR C30 is an electric suction unit for the aspiration of body liquids, oral, nasal and tracheal aspiration in adults or children. Provided with stand on 5 castors with brakes, it's suitable for hospital wards. Equipped with double jars, aspiration regulator and vacuum indicator located on the front panel. Footswitch control available upon request.

# ASKIR C30

Collection Jar in Autoclavable PC with Overflow Valve System

Antibacterial and Hydrophobic Filter

Vacuum Indicator (kPa and Bar)

Vacuum Regulator

Castors with Brake

Silicone Autoclavable Tubes



## TECHNICAL FEATURES

CE 0123

Motor	Oilless and maintenance-free piston pump
Power Feeding	230V-50Hz Other configurations available upon request
ISO 10079-1 Classification	HIGH VACUUM / HIGH FLOW
Max Vacuum (adjustable)	-0.80 Bar -80 kPa -600 mmHg
Max free air flow rate	40 l/min
Noise Level	61,5 dB
Power consumption	110 VA
Fuse	1 x F1.6 A 250V
Duty cycle	Non-stop operation
Weight	Kg 6.20
Size	cm 32 x 99 x 30
Years of Warranty	2
Shipping carton	1
Place of Manufacturing	Italy

## MAIN APPLICATIONS

Obstetrics

Endoscopy

General Surgery

Gastroenterology

Dental Practice

CA-MI

CA-MI srl Via Ugo La Malfa 31 - 43010 Pilastrò (PR) - ITALY  
Tel +39 0521 637133 - 631138 Fax +39 0521 639041  
export@ca-mi.it - vendite@ca-mi.it - [www.ca-mi.it](http://www.ca-mi.it)

SURGICAL  
SUCTION  
UNITS