



ECG-1215
ECG-1210 Digital 12-channel ECG

ECG-1215 Digital 12-channel ECG

- Features:**
- The foldable monitor(0~80 degree) with 7 inches color LCD;
 - 3×4+1, 6×2+1, 12, 12+1 waveforms display with grid line;
 - Real-time waveform freezing, replaying, printing;
 - Various analysis reports satisfy different requirement;
 - Optional SD card, unlimited data storage;
 - LAN port for real-time data transmission;
 - Multi-language options;
 - Long lasting Li-ion battery up to 2 hours continues working;
 - Easy upgrading keeps you updated;
 - Automatic arrhythmia delection and recording.



① Detailed Analysis Report



② LAN Port



③ Foldable Color LCD



④ ECG Information System



ECG-1210 Digital 12-channel ECG

- Features:**
- The foldable monitor(0~80 degree) with 5.7 inches LCD, contrast is adjustable;
 - R-R interval analysis, detailed and brief analysis reports satisfy different requirement;
 - Multi-language options;
 - Real-time waveform freezing;
 - Economy mode available help to reduce printing cost;
 - Long lasting Li-ion battery up to 2 hours continues working;
 - Optional SD card, unlimited data storage;
 - LAN port for real-time data transmission;
 - Easy upgrading keeps you updated;
 - Delightful operating experience with durable silicon button, unique shuttle design.



Specifications:

	ECG-1210	ECG-1215
Screen	5.7 inches LCD	7 inches color TFT
Acquisition Mode	Simultaneous 12-lead acquisition	
Input Impedance	>50 Ω	
Patient Leak Current	<10 μ A	
Time Constant	\geq 3.2s	
Frequency Response	0.05~150Hz	
Noise Level	<15 μ Vp-p	
Sensitivity	2.5, 5, 10, 20mm/mV	
Recording Paper Specification	210mm \times 140mm \times 20m fan-fold paper	
Paper Speed	5, 6, 25, 10, 12, 5, 25, 50mm/s	5, 10, 25, 50mm/s
Recording Format	3*4, 6*2+1, 6*2+3, 12, 12+1, Save mode, Arrhythmia mode	3*4, 3*4+3, 6*2, 6*2+1, 6*2+3, 12, 12+1, Cabrera mode, Save mode, Arrhythmia mode
Power Supply	AC: 100-240V, 50/60Hz, 60VA DC Sealed maintenance-free rechargeable battery continuous working time \geq 3 hours.	
Dimension	360mm \times 325mm \times 86mm	360mm \times 325mm \times 86mm
Weight	6.6Kg	8.6Kg



ECG-1215



ECG-1210