

**BECKER**



**MULTIPURPOSE MODULES FOR  
OR AND OUT-PATIENT DEPARTMENTS**

**INFORMATION**

**304  
E**

	DESCRIPTION	DIMENSIONS (W x D x H)	ORDER NO.	
	<p><b>Multipurpose modules for or and out-patient departments</b></p> <p>made of quality chrome nickel steel tubes and plates, alloy AISI 304 (1.4301), with polished surface to grain 240/320. All castors and legs conduct electricity.</p>			
	<p><b>Mayo instrument table</b> height manually adjustable 850 – 1280 mm, 4 swivel castors 75 mm Ø, 2 with brakes</p> <p>Tray size 600 x 400 mm</p> <p>Tray size 700 x 500 mm</p>		<p>130010</p> <p>130013</p>	
	<p><b>Oil-pump operated Mayo instrument table</b> oil hydraulically adjustable height 925 – 1310 mm, 3 castors 50 mm Ø</p> <p>Tray size 600 x 400 mm</p> <p>Tray size 700 x 500 mm</p> <p>Tray size 850 x 600 mm</p>		<p>130020</p> <p>130023</p> <p>130026</p>	
	<p><b>Oil-pump operated Mayo instrument table</b> oil hydraulically adjustable height 975 – 1360 mm, frame with removable tray, 1200 x 600 mm, 4 swivel castors 100 mm Ø, 2 with brakes</p> <p><b>Arm rest</b> height manually adjustable 790 – 1220 mm, without cushion, 5 swivel castors 75 mm Ø, 2 with brakes with fixed chrome nickel steel table top 325 x 450 mm</p> <p>with adjustable chrome nickel steel table top 325 x 450 mm</p> <p><b>Cushion</b> black, electrically conductive, for arm rest no. 130050 or 130055</p>	<p>1240 x 640 x 975 mm</p> <p>Ø 600 x 790 mm</p>	<p>130030</p> <p>130050</p> <p>130055</p> <p>130058</p>	
	<p><b>Hand and arm rest</b> height manually adjustable 680 – 980 mm, table top radiolucent, ball-type swivel joint tilts in all directions, without cushion, 4 swivel castors 75 mm Ø, 2 with brakes</p> <p><b>Cushion</b> black, electrically conductive, radiolucent, for hand and arm rest No. 130060</p>	<p>470 x 450 x 680 mm</p> <p>470 x 450 x 30 mm</p>	<p>130060</p> <p>130068</p>	

MULTIPURPOSE MODULES FOR  
OR AND OUT-PATIENT DEPARTMENTS

**BECKER**

	DESCRIPTION	DIMENSIONS (W x D x H)	ORDER NO.	
	<p><b>Kick bucket</b> without handle, 4 castors 50 mm Ø</p> <p>with chrome nickel steel bucket 15 l with plastic bucket, ivory, 10 l</p>	<p>450 x 450 x 425 mm</p>	<p>130110 130115</p>	
	<p><b>Kick bucket</b> with push handle, 4 swivel castors 50 mm Ø</p> <p>with chrome nickel steel bucket 15 l with plastic bucket, ivory, 10 l</p> <p><b>Chrome nickel steel lid</b> for chrome nickel steel bucket 15 l</p> <p>Chrome nickel steel bucket 15 l Plastic bucket, ivory, 10 l</p>	<p>450 x 450 x 825 mm</p>	<p>130120 130125  130127  130128 130129</p>	
 	<p><b>Single bowl stand</b> 1 chrome nickel steel bowl 8 l, push handle, with 4 castors 50 mm Ø with 5 swivel castors 75 mm Ø, 2 with brakes</p> <p>Spare bowl 8 l, chrome nickel steel</p>	<p>500 x 500 x 800/965 mm Ø 650 x 800/965 mm Ø 220/360 x 150 mm</p>	<p>130140 130141 130148</p>	
 	<p><b>Double bowl stand</b> 2 chrome nickel steel bowls 8 l, push handle, with 4 castors 50 mm Ø with 5 swivel castors 75 mm Ø, 2 with brakes</p> <p><b>Bowl stand with heating</b></p>	<p>500 x 790 x 800/965 mm 650 x 790 x 800/965 mm</p>	<p>130145 130146 130150</p>	
	<p><b>Linen hamper</b> with linen bag (ring 400 mm Ø) and push handle, with 4 castors 50 mm Ø with 5 swivel castors 75 mm Ø, 2 with brakes</p> <p>Spare linen bag</p> <p><b>Instrument receptacle</b> with push handle, without wire tray, with 4 castors 50 mm Ø with 5 swivel castors 75 mm Ø, 2 with brakes</p> <p><b>Perforated wire tray with folding handles</b></p>	<p>550 x 550 x 800/965 mm Ø 650 x 800/965 mm</p> <p>550 x 550 x 830/965 mm Ø 650 x 830/965 mm</p> <p>240 x 250 x 50 mm 480 x 250 x 50 mm</p>	<p>130160 130161 130168  130170 130171  120660 120665</p>	

	DESCRIPTION	DIMENSIONS (W x D x H)	ORDER NO.	
	<p><b>Surgeon's step</b> removable, grooved non-slip plate</p>	<p>600 x 350 x 130 mm 600 x 350 x 220 mm 600 x 350 x 300 mm</p>	<p>130200 130203 130206</p>	
	<p><b>Surgeon's double step</b> removable, grooved non-slip plates</p>	<p>600 x 650 x 220/440 mm</p>	<p>130210</p>	
	<p><b>Wall-mounted system to store 4 surgeon steps</b> Pair of rails for wall mounting with hooks to take max. 4 surgeon steps no. 130200, 130203 and/or 130206; without surgeon steps</p>	<p>600 x 430 x 1450 mm</p>	<p>130220</p>	
	<p><b>Infusion stands</b> height adjustable 1370 – 2500 mm, 4 hooks</p> <p>with PP drip catcher, 5 swivel castors 75 mm Ø, 2 with brakes</p> <p>with PP drip catcher, 5 castors 50 mm Ø</p> <p>„Model Ulm“ without drip catcher, with push handle, 5 swivel castors 75 mm Ø, 2 with brakes</p>	<p>650 x 650 x 1370 mm</p>	<p>130240 130243 130246</p>	

MULTIPURPOSE MODULES FOR  
OR AND OUT-PATIENT DEPARTMENTS

**BECKER**

	DESCRIPTION	DIMENSIONS (W x D x H)	ORDER NO.
	<p><b>Accessory rack</b> 2 standard rails 25 x 10 mm and 1 basket 1/2-STE, 4 swivel castors 75 mm Ø, 2 with brakes</p>	600 x 490 x 1450 mm	130260
	<p><b>Infusion preparation table</b> with six infusion hooks, fixed chrome nickel steel work top, 2 baskets 1/2-STE with opening, 4 swivel castors 75 mm Ø, 2 with brakes</p>	600 x 500 x 1600 mm	130270
	<p><b>Infusion preparation table</b> with ten infusion hooks, fixed chrome nickel steel work top, 2 storage baskets 900 x 265 mm, 4 swivel castors 75 mm Ø, 2 with brakes</p>	910 x 600 x 1600 mm	130275
	<p><b>Anaesthesia accessory rack</b> with six hooks, panel with 6 metal rods and special tray, 4 swivel castors 75 mm Ø, 2 with brakes</p>	600 x 500 x 1600 mm	130280
	<p><b>Anaesthesia accessory transport trolley</b> 3 Pairs of hooks for respiration tubes, 2 hooks for respiration bags, without tubes and bags, guard frame on the top, at the back 4 pairs of hooks with 2 baskets 1/2-STE and 2 baskets 1-STE with opening, 4 swivel castors 75 mm Ø, 2 with brakes</p>	600 x 680 x 1920 mm	130285

	DESCRIPTION	DIMENSIONS (W x D x H)	ORDER NO.	
	<p><b>Transport trolley</b> for linen or waste bags, with push handle, without lid; bottom plate, removable, 4 swivel castors 75 mm Ø, 2 with brakes</p> <p>single</p> <p>double</p> <p><b>Chrome nickel steel lid</b> for transport trolley No. 130300 or 130305</p>	<p>465 x 490 x 900 mm</p> <p>770 x 490 x 900 mm</p>	<p>130300</p> <p>130305</p> <p>130308</p>	
	<p><b>Operating room shoe rack</b> with 6 hanging wire baskets for 24 pairs of shoes</p> <p>mobile on 4 swivel castors 75 mm Ø, 2 with brakes</p> <p>for wall-mounting</p>	<p>910 x 600 x 1600 mm</p> <p>910 x 280 x 1400 mm</p>	<p>130310</p> <p>130315</p>	
	<p><b>Paper rack</b> for 4 different types of paper up to 1300 mm length, 4 swivel castors 100 mm Ø, 2 with brakes</p>	<p>1400 x 650 x 920 mm</p>	<p>130320</p>	
	<p><b>Chrome nickel steel shelf trolley</b> with fix welded shelves, 4 swivel castors 125 mm Ø, 2 with brakes</p> <p>with 3 shelves</p> <p>with 4 shelves</p>	<p>1340 x 670 x 1170 mm</p> <p>1340 x 670 x 1600 mm</p>	<p>130330</p> <p>130335</p>	

MULTIPURPOSE MODULES FOR  
OR AND OUT-PATIENT DEPARTMENTS

**BECKER**

	DESCRIPTION	DIMENSIONS (W x D x H)	ORDER NO.	
	<p><b>Container transport trolley</b> for sterile material containers 600/300 mm, 4 swivel castors 125 mm Ø, 2 with brakes</p> <p>unilateral with 6 chrome nickel steel trays</p> <p>bilateral with 12 chrome nickel steel trays</p>	<p>650 x 450 x 1540 mm</p> <p>650 x 750 x 1540 mm</p>	<p>130350</p> <p>130355</p>	
	<p><b>Hook bar for OR aprons</b> for wall-mounting double hooks, movable</p> <p>500 mm with 3 double hooks</p> <p>600 mm with 4 double hooks</p> <p>800 mm with 5 double hooks</p> <p>1000 mm with 6 double hooks</p> <p>1200 mm with 8 double hooks</p> <p>Spare double hooks</p>		<p>130360</p> <p>130362</p> <p>130364</p> <p>130366</p> <p>130368</p> <p>130378</p>	
	<p><b>Bed screen mobile</b> 3 swivel castors 75 mm Ø with brakes, without curtain</p> <p><b>Bed screen fixed</b> without curtain</p> <p><b>Plastic curtain</b> white</p> <p><b>Linen curtain</b> white</p>	<p>900 x 1900 mm</p> <p>900 x 1800 mm</p>	<p>130390</p> <p>130395</p> <p>130397</p> <p>130398</p>	

# BECKER



**Becker quality steel furniture and equipment have been in use in hospitals and medical centers for decades. Optimum use of available space is achieved through an extensive line of programmes and sizes.**

**Special designs and individual choice of colors as well as custom-built systems offer an enormous variety of options to the customer.**

**Metallwarenfabrik  
Walter H. Becker GmbH  
Anzenkirchener Straße 4  
D-84371 Triftern**

**Tel.: +49 (0) 85 62 / 96 05-0  
Fax: +49 (0) 85 62 / 24 34  
E-mail: becker.triftern@t-online.de  
Internet: www.becker-triftern.de**

**BECKER SYSTEMS for:**

**Wards  
Out-patient Departments  
Operating Rooms  
Treatment Rooms  
Sterilization**