



DCP-900

Auto Chart Projector



DONGYANG Optics Co., Ltd.

DCP-900

Auto Chart Projector

THE NEW CHART PROJECTOR DCP-900
WITH VARIOUS CHARTS OFFERS INNOVATIVE
FUNCTIONS FASTER THAN EVER BEFORE.


DONG YANG's New Auto Chart Projector DCP-900 is our successful achievement for refraction examination with maximum care for high resolution image.

New Auto Chart Projector DCP-900 provides you with the function of which remote controller is available to select one of the signal channels in order to prevent other signals from receiving.


Specifications	
Chart	31pcs
Chart selection speed	Average 0.3 sec
Mask	1 open / 3 horizontal lines / 3 vertical lines / 21 single letters
Filter	Red/Green
Projection distance	2.0m~7.0m(3m is standard)
Projection magnification	30 x (at 5m)
Projection size	ø 252mm and 330W x 225H (at 5m)
Tilt angle	Ball joint ($\pm 10^\circ$)
Lamp	6V 20W halogen
Power source	AC100~240V 50/60Hz
Power consumption	80VA
Power saving function	5~10minutes after
Dimensions	Main body 265(H) x 230(W) x 360(D)mm Remote controller 20(H) x 64(W) x 195(D)mm
Net weight	Main body 6.0Kg / Remote controller 180g
Interface	8-pin connector for simultaneous control of chart, mask and R/G filter by a motorized refractor's control unit.
Accessories	1 each of dust cover 1, Polarized metal screen(400mm x 350mm), 1 polarized glasses, 2 spare fuses and 2 batteries, 1 spare lamp
Optional accessories	Wall bracket, red/green glasses

Optotypes



A type



B type



C type



DONGYANG Optics Co., Ltd.

#1101-1, 6th Daeryung Techno Town, Gasan-Dong, Geumcheon-Gu, Seoul, Korea
Tel +82-2-7573234 Fax +82-2-754-8189 www.dyoapt.com E-MAIL : trade@dyopt.com