




- 01 | ARTHROSCOPY
- 02 | GYNECOLOGY
- 03 | LAPAROSCOPY
- 04 | ENT
- 05 | THORACOSCOPY**
- 06 | UROLOGY
- 07 | BRONCHOSCOPY
- 08 | PAEDIATRY
- 09 | SCOPE WITH
ANGLED EYEPIECE
- 10 | ACCESSORIES

THORACOSCOPY

Lösungen für das medizinische Untersuchungs- und
 Therapieverfahren der Thorakoskopie

Solutions for the medical examination and
 treatment procedures of thoracoscopy



Scale: 1:2
 Maßstab: 1:2

Article number Artikelnummer	Ø Diameter Ø Durchmesser	Direction of View Blickrichtung	Working Length Nutzlänge
EB-0290.00.360	2,9 mm	0°	360 mm
EB-0290.30.360	2,9 mm	30°	360 mm
EB-0650.30.360	6,5 mm	30°	360 mm
ET-1000.00.175	10,0 mm	0°	175 mm
ET-1000.00.225	10,0 mm	0°	225 mm
ET-1000.30.175	10,0 mm	30°	175 mm
ET-1000.30.225	10,0 mm	30°	225 mm

Available Directions of View
 Verfügbare Blickrichtungen

0°

30°

