

PLUG & TOUCH BY RODENSTOCK

PD MAX
DIGITAL PUPILOMETER

Digital. Direct-acting. Ergonomic.



- Internal and external readings
- Energy saving by auto power off
- Illuminated internal LCD display
- Variable measuring distance
- Measuring principle: cornea light reflexion

QUALITY IN DETAIL



Simultaneous data display

The compact **Rodenstock PD Max** digital pupilometer has the unique advantage of having two separate LCD displays. In addition to the large external LCD display the **PD Max** also has an integrated illuminated LCD. Within the unit, in line with the viewing field, the user can simultaneously see the results during the measurement.

Precision

The precise "cornea light reflection" measurement system provides fast and reliable results with a convenient adjustable focal distance from 30 cm to infinity.

High usability

The design of the **PD MAX** offers the user an intuitive and easy to use tool. Not only the usability but also the selection of the internal light level and auto off function make this device to a perfect solution.

SPECIFICATIONS

Measurement

Binocular	46.0 to 83.0 mm
Monocular	23.0 to 41.5 mm
Distance	30 mm to infinity
Steps	30, 35, 40, 50, 62 cm, 1.0, 2.0 m, infinity
Appreciation	0.5 mm
Precision	1.0 mm

Display

Fixation

Fixation lamp

Occlusion

Dimensions WDH

Standard

Battery type



2012/04 - subject to change without notice