


ComfortStar® - Anesthesia face masks Consumables and Accessories

Patient comfort and ease of use were the top priorities when developing the Dräger disposable ComfortStar® anesthesia face masks. Patient comfort is supported by a teardrop or a round-shaped mask, which provides an anatomical fit with effective seal. Non-slip ridges and a flexible mask body provide healthcare professionals an easy and comfortable grip on the mask.



Benefits

ComfortStar anesthesia face masks are for disposable use

The inflatable mask cushion and the flexibility of the mask give the physician a comfortable grip. In addition, the ComfortStar masks sizes 4-6 are equipped with slip-stop ridges to provide a safe and non-slip grip.

The transparent mask body allows for continuous observation of the patient's condition. All masks are pre-equipped with 'Traffic Light' colored hook rings. They are color coded for easier size identification. Smaller sizes (0-1) have reddish color, mid-sizes (2-3) yellowish color and large sizes (4-6) greenish color. Small and large headgears are also available for easy and secure fixation.

A potential side effect for patients who have undergone surgical procedure under general anesthetic is post-operative nausea and vomiting (PONV). This often delays patient recovery. Research shows that added aromas can help minimizing these unpleasant side effects¹⁾.

ComfortStar scented anesthesia face masks have the same excellent benefits of the ComfortStar standard masks, but in addition have added aromas. They are available with strawberry or bubble gum scents in sizes 1-4, and mint scent in sizes 3-6. The appropriate color scheme provides optical stimulation and is intended to reduce pre-operative anxiety in patients.

¹⁾ Anderson, L.A., Gross, J.B.: Aromatherapy with peppermint, isopropyl alcohol, or placebo is equally effective in relieving postoperative nausea, 2004, in: Journal of PeriAnesthesia Nursing, Vol. 39, No 1, 29-55

Compatibility

All ComfortStar masks are compatible with all Dräger anesthesia devices.

Ordering Information

Description	Package of	Order No.
ComfortStar anesthesia face mask, disposable, standard, with hook ring		
size 0, neonate	20 pcs.	MP01510
size 1, new born	20 pcs.	MP01511
size 2, infant	20 pcs.	MP01512
size 3, toddler	20 pcs.	MP01513
size 4, small adult	20 pcs.	MP01514
size 5, medium adult	20 pcs.	MP01515
size 6, large adult	20 pcs.	MP01516
ComfortStar anesthesia face mask, disposable, strawberry scented, with hook ring		
size 1, new born	20 pcs.	MP01521
size 2, infant	20 pcs.	MP01522
size 3, toddler	20 pcs.	MP01523
size 4, small adult	20 pcs.	MP01524
ComfortStar anesthesia face mask, disposable, bubble gum scented, with hook ring		
size 1, new born	20 pcs.	MP01531
size 2, infant	20 pcs.	MP01532
size 3, toddler	20 pcs.	MP01533
size 4, small adult	20 pcs.	MP01534
ComfortStar anesthesia face mask, disposable, mint scented, with hook ring		
size 3, toddler	20 pcs.	MP01543
size 4, small adult	20 pcs.	MP01544
size 5, medium adult	20 pcs.	MP01545
size 6, large adult	20 pcs.	MP01546
Headgear, Silicone, reusable		
small	1 pc.	MP01557
large	1 pc.	MP01558
Headgear, Soft Fabric, disposable		
small	1 pc.	MP01559
large	1 pc.	MP01560
Headgear, Soft Rubber, disposable		
small	1 pc.	MP01561
large	1 pc.	MP01562
Two-sided pump ball for air management	1 pc.	MP01590

Notes

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23558 Lübeck, Germany
www.draeger.com

Manufacturer:

Dräger Medical GmbH
Moislinger Allee 53-55
23558 Lübeck, Germany

As of August 2015

Dräger Medical GmbH changes
to Drägerwerk AG & Co. KGaA

Locate your Regional
Sales Representative at:
www.draeger.com/contact



REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH
Moislinger Allee 53-55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION ASIA / PACIFIC

Draeger Medical
South East Asia Pte Ltd.
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION MIDDLE EAST, AFRICA

Dräger Medical GmbH
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama Comercial
S. de R.L.
Complejo Business Park,
V tower, 10th floor
Panama City
Tel +507 377 9100
Fax +507 377 9130
contactcsa@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute Technologie
d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
d1mfr-contact@draeger.com

REGION NORTH AMERICA

Draeger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free+1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

The quality management system at Dräger Medical GmbH is certified according to ISO 13485, ISO 9001 and Annex II.3 of Directive 93/42/EEC (Medical devices).