

PRIME VERSIONS WITHOUT AIR DRYER









50/15 PRIME S







Code: PR-0F200-040

Description: Oil free compressors ,FAD at 5 bar 170 l/min, tank capacity 40 l., without air dryer.

This compressor is available in 10 BAR VERSION also! visit the "10 bar oil-free compressors section" for more details!

Model									
	Dental chairs	Air dryer	Air Capacity m ³ /h l/min cfm	FAD @ 5bar m ³ /h l/min cfm	Motor HP KW	Volt. V	Freq. Hz	Receiver l GAL	
50/15 PRIME S	2-3	no	- 255 9,01	- 170 6	2 1,5	230	50	40 10,57	

Model		Duty cycle			
	Noise Level dB(A)		Max pressure bar PSI	Packaging mm	Net weight Kg lb
50/15 PRIME S	66	70%	7 101	600x410x770	46 101

