

PRIME VERSIONS WITHOUT AIR DRYER








50/25 PRIME S







Code: PR-OF300-040

Description: Oil free compressors, 3 cylinders, FAD at 5 bar 250 l/min, tank capacity 40 l, without air dryer, for heavy duty cycle.

This compressor is available in 10 BAR VERSION also! visit the "10 bar oil-free compressors section" for more details!

Model							
	Dental chairs	Air dryer	Air Capacity m ³ /h l/min cfm	FAD @ 5bar m ³ /h l/min cfm	Motor HP KW	Volt. Freq. V Hz	Receiver l GAL
50/25 PRIME S	4	no	- 360 12,7	- 250 8,8	3 2,2	230-400 50	40 10,57

Model		Duty cycle			
	Noise Level dB(A)		Max pressure bar PSI	Size mm	Gross weight Kg lb
50/25 PRIME S	66	70%	7 101	600x410x810	54 119

