

PRIME VERSIONS WITHOUT AIR DRYER

270/100 TANDEM PRIME S














Code: PR-OF1200-270-T3

Description: The new models 270/100 PRIME are equipped with the new 6 cylinders PRIME motors.

This motor is the result of the constant MGF research and development activity, able to provide 900 l/min at 5 bar.

This compressor is available in 10 BAR VERSION also! visit the "10 bar oil-free compressors section" for more details!

Model							
	Dental chairs	Air dryer	Air Capacity	FAD @ 5bar	Motor	Volt. Freq.	Receiver
			<i>m3/h l/min cfm</i>	<i>m3/h l/min cfm</i>	HP KW	V Hz	<i>l GAL</i>
270/100 TANDEM PRIME S	15	no	- 1.440 50,8	900 31,6	15 10	400-50 50	270 71.3

Model				
	Noise Level	Duty cycle	Max pressure	Size
	<i>dB(A)</i>		<i>bar PSI</i>	<i>mm</i>
270/100 TANDEM PRIME S	75	70%	7 101	1560x1000x1030
				<i>Net weight</i>
				<i>Kg lb</i>
				220 485

