

PRIME VERSIONS WITH AIR DRYER



270/100 TANDEM PRIME M







Code: PR-OF1200-270-MT3

Description: The new models 270/100 PRIME are equipped with the new 6 cylinders PRIME motors and membrane dryer system, able to provide a dew point value till -20°C. The model 270/100 PRIME provides 780 l/min at 5 bar.

This compressor is available in 10 BAR VERSION also! visit the "10 bar oil-free compressors section" for more details!

Model							
	Dental chairs	Air dryer	Air Capacity	FAD @ 5bar	Motor	Volt. Freq.	Receiver
			m ³ /h l/min cfm	m ³ /h l/min cfm	HP KW	V Hz	l GAL
270/100 TANDEM PRIME M	14	yes	- 1440 50,8	780 27,5	13 10	400-50 50-60	270 71,3

Model				
	Noise Level	Duty cycle	Max pressure	Size
	dB(A)		bar PSI	mm
270/100 TANDEM PRIME M	75	70%	7 101	1560x1000x1030
				Net weight
				Kg lb
				240 529

