









GENESI versions with air dryer






24/10 GENESI M



Code: D-GE150-024-M

Description: Oil free compressor, FAD at 5 bar 120 l/min, tank capacity 24 l., with air dryer .

Model	 Dental chairs	 Air dryer	 Air Capacity	 FAD @ 5bar	 Motor	 Volt.	 Freq.	 Receiver
			<i>m3/h l/min cfm</i>	<i>m3/h l/min cfm</i>	<i>HP KW</i>	<i>V</i>	<i>Hz</i>	<i>l GAL</i>
24/10 GENESI M	1-2	yes	- 182 4,2	- 120 4,2	1,5 1,1	230	50	24 6.34

Model	 Noise Level	 Duty cycle	 Max pressure	 Size	 Net weight
	<i>dB(A)</i>		<i>bar PSI</i>	<i>mm</i>	<i>Kg lb</i>
24/10 GENESI M	65	50%	7 101	420x500x600	33 72

