

Silent PRIME versions with air dryer




CS 100/50 TANDEM PRIME M







Code: PR-CS600-090-MT3

Description: Silenced oil free compressor, 6 cylinders, FAD at 5 bar 450 l/min, tank capacity 90 l, with membrane air dryer able to warranty till -20° C dew point value, compressor suitable for heavy duty cycles

This compressor is available in 10 BAR VERSION also! visit the "10 bar oil-free compressors section" for more details!

Model							
	Dental chairs	Air dryer	Air Capacity	FAD @ 5bar	Motor	Volt. Freq.	Receiver
			m ³ /h l/min cfm	m ³ /h l/min cfm	HP KW	V Hz	l GAL
CS 100/50 TANDEM PRIME M	7	yes	- 720 25,4	- 450 15,8	3+3 2,2+2,2	230-400 50	90 23,77

Model				
	Noise Level	Duty cycle	Max pressure	Size
	dB(A)		bar PSI	mm
CS 100/50 TANDEM PRIME M	56	70%	7 101	725x1245x1020
				Net weight
				Kg lb
				225 495

