

# medel®

## CHECK

### BLOOD PRESSURE MONITOR



Universal cuff arm  
22-42 cm



Extra-large backlight display



Irregular heartbeats indicator



Average of the 3 measured



2 users



5 year warranty

Medel **CHECK** is a fully automatic arm blood pressure monitor; it is supplied with T3 function that allows to read the average of the last 3 measurements taken in 15 minutes. Medel CHECK comes with a large and backlight display for an easy reading. Tested in accordance with ESH (European Society of Hypertension) Medical Device Class IIa.

#### T3 function

it simply and automatically displays the average of the last 3 readings

#### Ergonomic handle



#### 2 Users

with a memory of 60 measurements for each user

# CHECK

Medel CHECK is equipped with a universal cuff. It is supplied with a sturdy case for travelling and transportation.



Technical specifications	CHECK
Technology	Oscillometric method
Accuracy	Pressure +/- 3 mmHg Heart rate +/- 5%
Measurement range	Pressure 30 to 280 mmHg; heart rate 40 to 199 beats/min.
Power supply	4 AA 1,5V batteries and USB port
Memory	2 users: 60 memories for each user
Systolic/Diastolic and Pulse display	√
Date and time	√
Irregular heartbeat indicator (IHB)	√
T3 Function	√
LCD size	6 x 8 cm
Cuff Type	Universal 22-42 cm

COD. 92821 Medel CHECK

CE 0123

