

medel®

iCARE

Automatic blood pressure monitor with *Bluetooth*®



60 Memory spaces for 2 Users



Medel Health Check app



Data transmission on smartphone



Classification of the measured values



5-Year Warranty

Medel **iCARE** is the upper arm blood pressure monitor used to carry out the measurement and monitoring of the arterial blood pressure values. This allows to quickly and easily measure your blood pressure values. Data can be transferred on your smarhnone using the " Medel Health Check" app, available on Appstore and Google Play

- 2 Users x 60 memory spaces • Average of the measured values (morning and evening) • Data transfer by Bluetooth® technology • Classification of values with coloured scale • Arrythmia alert • Date and time of measurements • Storage pouch and batteries included




Universal cuff
22-44 cm



Date and time of
measurements



Extra-large display



TECHNICAL SPECIFICATIONS

Technology	Oscillometric
Accuracy	Pressure +/- 3 mmHg, Heart beats +/- 5%
Measurement range	Pressure from 30 to 250 mmHg; Pulse 40/180 heart beats/min.
Power supply	4 1,5 V AAA alkaline batteries
Memory	2x60 posizioni di memoria
Users	2
Systolic/Diastolic and Pulse display	✓
Date and time	✓
Irregular heart beat indicator	✓
Average of measurements	✓
LCD size (bxh)	7,6 x 7 cm
Cuff type	Soft
Cuff size	Universal 22-44 cm



ACCESSORIES

- Universal cuff 22-44 cm
- Storage pouch

95164

