

4500 MRI Pulse Oximeter

Fiberoptic SpO₂ Monitoring for MRI



Operational Features:

- Fiberoptic SpO₂ Technology
- Heart Rate and SpO₂ Values
- Adult, Pediatric & Neonatal
- Grip Sensor Compatible
- Bright LED "Smart Display"
- Battery Operation

Specifications:

Accuracy

100-90%, ± 2.0%
90-80%, ± 2.2%
80-70%, ± 2.5%

Features

Pulse Rate Range 20 - 250 bpm

Product Dimensions & Weight

Height: 3.25 in (8.25 cm)
Width: 6.63 in (16.8 cm)
Depth: 10.25 in (26.03 cm)
Weight: 5.6 lbs (3.15 kg)

Provides continuous and accurate
monitoring for increased patient safety
during MRI procedures



12601 Research Parkway
Orlando, FL 32826
tel +1-407-275-3220
tel 1-800-331-3220
fax +1-407-249-2022
info@invivocorp.com

www.invivocorp.com

MADE IN THE U.S.A.

© 2006 Invivo Corporation All Rights Reserved
580099 Rev 01