



Discreet. Dignified. Dependable.

2015 Product Catalog







Table of Contents

Acute Care

Cath Grip® - Catheter Securement System	
UniGrip™	3
Grip & Seal™	4
Hold Me™	5
Safe n' Dry®	6

Skin Interface Products

BioPlus+™	7
FreeDerm™	7

Urine Drainage Bags

Urine Drainage Leg Bag 1,000mL	8
Urine Drainage Bedside Bag 2,000mL	8
Urine Drainage Bedside Bag 2,000mL with Urine Meter	8

Direct to Consumer

Men's Liberty	9
----------------------------	---

International Products

XLS External Continence Device for Men™	10
Cath Grip® UniGrip™ -	10

Product Summary	11
------------------------------	----

Contact Information	12
----------------------------------	----





CathGrip provides universal securement for catheters from 2-47 Fr. With a gentle, skin-friendly hydrocolloid base, CathGrip can help eliminate skin tears, and moves with the body for improved skin outcomes and comfortable patients.

UniGrip™


UniGrip is a truly universal system that securely holds tubes from 2-47 Fr. Available in single or dual straps for all over securement with a skin friendly hydrocolloid adhesive to prevent skin injuries.

Soft straps with no-slip grip, will work for both adult and pediatrics. Gentle grips are durable for repeat positioning.


Available as sterile or non-sterile.



Uses listed above are not all inclusive. May vary based on tube location, infection prevention protocols and presence of exuding wounds.



***Each SKU fits
2-47 French!***





CathGrip provides universal securement for catheters from 2-47 Fr. With a gentle, skin-friendly hydrocolloid base, CathGrip can help eliminate skin tears, and moves with the body for improved skin outcomes and comfortable patients.


Grip & Seal™

Grip & Seal extends the security of our hydrocolloid base to support insertion site securement.

Hydrocolloid wings form an occlusive seal that securely holds any tube from 2-47 Fr and can help reduce dislodgement and infection rates.

Only available as sterile.


Suggested Grip & Seal Applications



Grip & Seal™ is only recommended for well established tubes that do not require regular rotation. Do NOT use on exuding wounds.



Insertion Site Securement for Tubes 2-47 Fr





CathGrip provides universal securement for catheters from 2-47 Fr. With a gentle, skin-friendly hydrocolloid base, CathGrip can help eliminate skin tears, and moves with the body for improved skin outcomes and comfortable patients.



HoldMe™

Designed to prevent pressure injuries and skin tearing with low-profile proprietary hydrocolloid securement. Especially important when dealing with tender skin seen in senior or pediatric populations.

Low profile oval pairs that are sized for difficult to apply areas. Will hold securely for up to 7 days.

Available as sterile or non-sterile.

Suggested Hold Me Applications



Low Profile Securement for Hard to Reach Areas



Safe n' Dry™

The completely external Safe n' Dry seals gently and securely to the tip of the penis with BioDerm's proprietary hydrocolloid adhesive XtraDerm™, designed for institutional and hospital use.

Safe n' Dry:



- Has a completely external design
- Fits most male anatomy
- Keeps men dry and comfortable
- Allows men to intermittent cath with the device in place
- Average wear time of 24-48 hours
- Replaces diapers, condom catheters and medically unnecessary Foley catheters



Clinical protocol available for institutional use.

Ideal for:

- Hospitals
- Rehab Center
- Long Term Acute Care
- Nursing Homes
- Other In-Patient Facilities



Skin Interface Products



FREEDERMTM Adhesive Remover

Instantly removes hydrocolloid without irritation, evaporates rapidly with no residue.

Allows immediate re-application of an adhesive device without washing. Available as 1 oz pump spray or boxes of 25 wipes.



Bio Plus+ Barrier Prep Wipe

Improves hydrocolloid adhesion by up to 50%. A BioPlus+ wipe is included in most BioDerm products and is also available as a stand alone product.



Drainage Bags

1000ml Leg Bag

Features:

- Holds up to 1000mL
- Kink-resistant tubing
- Convenient drain spout
- Anti-reflux valve
- All fluid pathways are sterile
- Leg bag straps included



2000ml Bed Bag

Features:

- Holds up to 2000mL
- Kink-resistant tubing
- Convenient drain spout
- Anti-reflux valve
- All fluid pathways are sterile
- Bedside hooks included



2000ml Bed Bag w/Urine Meter

Features:

- Holds up to 2000mL, plus 100mL urine meter for monitoring
- Kink-resistant tubing
- Convenient drain spout
- Anti-reflux valve
- All fluid pathways are sterile
- Bedside hooks included





The First New Technology for Managing Men's Incontinence in 50 Years!

The completely external Men's Liberty seals gently and securely to the tip of the penis with BioDerm's proprietary hydrocolloid adhesive XtraDerm™, designed for extended wear time. Men's Liberty keeps the skin dry and comfortable—all day, every day.



How Does It Work?

Men's Liberty is made from hydrocolloid, a medical adhesive that creates a secure, skin friendly seal around the urinary opening. Hydrocolloid has been used successfully in wound and ostomy treatment for decades.

Men's Liberty:

- Has a completely external design
- Fits most male anatomy
- Keeps men dry and comfortable
- Allows men to intermittent cath with the device in place
- Average wear time of 24-48 hours
- Replaces diapers, condom catheters and medically unnecessary Foley catheters
- Covered by Medicare, most state Medicaid, Private Insurances, Workers Compensation and VA/TriCare

HCPCS Code: A4356



External urethral compression device. One size fits all. Designed to temporarily halt the flow of urine through the urethra while applying Men's Liberty™.

HCPCS Code: A4356




International Models

XLS[®]

BioDerm's XLS provides the same features as Men's Liberty. It does not have the integral collection chamber and connects directly to a leg or bedside bag.

Only for sale outside US.

Currently available in: Canada, Great Britain, Ireland, Germany, the Netherlands and Peru.




CathGrip[®]

BioDerm's international universal securement option:

- Retains all the key features of the UniGrip securement system
- Only available as non-sterile

Provides safe, simple securement for:

- Urinary
- PEGs
- Drains
- Chest Tubes
- And many more...





Product Summary

Product Name	HCP Code	Catalog/SKU Number*	Size (if applicable)	Units Per Box	Boxes Per Case
UniGrip™					
UniGrip Large	A4333 for Urinary	51300	4"x4.75"	30	3
UniGrip Large (1 Strap)		51350	4"x4.75"	30	3
UniGrip Medium		51301	3"x3.5"	30	3
UniGrip Small		51250	3"x3.5"	30	3
Grip & Seal™					
Grip & Seal Large	A5200 for Percutaneous	51305	4.32"x7"	30	3
Grip & Seal Medium		51303	2.94"x4.75"	30	3
Grip & Seal Small		51100	2.2"x3.56"	30	3
HoldMe™					
HoldMe X-Large	Varies Based on Type of Tubing Secured	54009	4.32"x7"	50	6
HoldMe Large		54003	3"x3.5"	25	12
HoldMe Medium		53002	1.48"x2.4"	25	12
HoldMe Small		54002	.93"x1.5"	25	12
HoldMe X-Small		54001	.62"x1"	25	12
Skin Interface Products					
BioPlus+™ Wipe	A5120	52300	One Size	25	12
FreeDerm™ Spray	A4455	52200	1 oz.	8	4
FreeDerm™ Wipe	A4456	52201	One Size	25	8
Urological					
Men's Liberty™	A4326	23046	One Size	35	12
Penile Clamp™	A4356	20003	One Size	1	10
Safe n' Dry™	NA	21026	One Size	30	12
International Models					
XLS™	NA	21026	One Size	10	36
UniGrip™ (non-sterile only)	NA	51400	One Size	30	3
Drainage Bags					
Leg Bag	A4358	90000	1000mL	50	4
Bedside Bag	A4357	90001	2000mL	50	4
Bedside Bag with Urine Meter	NA	90005	2000mL	10	4

*Add suffix NS to UniGrip or HoldMe to order non-sterile product.

For more information please visit us online at www.BioDermInc.com



Contact Us

If you have any questions, would like more information or to receive a sample of one of our products, please contact our Customer Care Team or visit us online to submit a request form.

BioDerm, Inc.

12320 73rd Court North
Largo, FL 33773

Customer Service:

(800) 814-3174

CustomerCare@BioDermInc.com

Check Out The Men's Liberty Blog!

We are publishing a new incontinence blog with excellent information, news and research on incontinence, related diagnoses and caregiving. We welcome guest blogs from physicians and other health care professionals.

Visit Us Online for More Information at:

www.BioDermInc.com

Find us on Social Media!

