


Kodak | Dental Systems

Kodak CR 7400
Digital Radiography System

Kodak CR 7400

Digital Radiography System



Kodak CR 7400

Digital Radiography System

A Smooth Digital Transition

The Kodak CR 7400 digital radiography system is an easy way to transition your dental practice to digital radiography using your current film radiography equipment, in most cases – including your intraoral, panoramic and cephalometric devices. It's a quick and simple transition to creating digital intraoral and extraoral radiographic images that can be viewed and stored on your computer.



An Easy, Familiar Workflow

Of all the digital radiography technologies, a “computed radiography” (or CR) system has a workflow most like film. Not only is it easy to learn how to use, it also causes minimal disruption in your workflow and your office layout.

Expose an image using the same x-ray equipment, positioners and techniques that you use for film radiography. Instead of film, you use a radiation-sensitive CR imaging plate.

1

Kodak CR 7400
Digital Radiography System

Intuitive software interface allows you to quickly view images on a computer monitor, either for diagnosis or for sharing with your patient. Images are then stored on your computer.

3

Digital made simple.

Scan from one image to a series of images, including FMS, in one step using the compact table-top scanner. Plates are automatically erased to save time and prevent double exposures.

2



Larger, high quality images



An Affordable, Cost-effective Digital Option

While many digital radiography systems may require purchasing new equipment, the Kodak CR 7400 system will most likely work with your current film x-ray equipment. As a result, the initial investment in dollars, time and training for the transition to a CR system is much less than other digital radiography options. It's also easy to share one Kodak CR 7400 system between multiple operatories from a centralized location.

Improved Patient Care

Most importantly, the Kodak CR 7400 system gives you the improved diagnostics and patient care opportunities of digital radiography. The high quality images can be enlarged far beyond the size of film for easy viewing. And the Kodak dental imaging software allows you to enhance and manipulate images for diagnostic advantages.


Many dental professionals show on-screen images to patients in the operator to educate and explain treatment procedures. And since the plates are flexible and wireless, and come in sizes 0 to 4, many patients find them more comfortable than the thicker, rigid sensors of other direct digital systems.

The Flexible Digital Practice

Not only can you share one Kodak CR 7400 system between operatories, but you can add other Kodak digital systems to your practice, confident that they will seamlessly integrate with Kodak imaging and practice management software. Many practices desire a mix of computed radiography (CR) and direct digital radiography (DR) systems to take advantage of their different benefits.



One-step scanning



Durable and
dependable

Proven, Dependable Kodak Technology

Behind every Kodak CR 7400 system is Kodak's reputation for dependability. Our imaging plates are durable, and replacement plates are relatively inexpensive. Plates are fitted with a light protection Kodak pouch, which guards against light damage and scratches. For each use, plates are slipped into a disposable hygienic barrier envelope. Kodak's thin and flexible imaging plates offer low noise and grain for high image quality.

Take the Next Step

Because Kodak offers the entire range of radiographic solutions, from digital to film, we can help you figure out exactly which systems are right for your practice.

If you'd like more information, or if you'd like to schedule a demonstration, call our sales and service center. (Please refer to the back page for our free telephone number)

To learn more about Kodak's complete suite of digital imaging systems, go to www.kodak.com/dental

Support and Guidance You Can Count On

Kodak is committed to making your transition to digital radiography as smooth and seamless as possible, including after your purchase.

The Kodak CR 7400 system comes with:

- No ongoing routine equipment maintenance
- Optional replacement program to minimize any downtime
- Unlimited toll-free telephone and email support*

*Warranty and service plans vary by country.



Would you like to know more?

To schedule a demonstration or to receive further information,
please call 0800 89 51 13 (United Kingdom, toll-free number)
1800 55 49 97 (Ireland) - 800 1 71 49 (Denmark),
0800 11 49 01 (Finland) - 800 1 11 67 (Norway)
020 79 74 25 (Sweden) - +49 (0)711 406 5674 (Other countries)
or visit our website: www.kodak.com/dental



© Eastman Kodak Company, 2006.
KODAK is a trademark of Eastman
Kodak Company.