



- Complete protection against dust and liquids IP68
- Washable, sterilized in just seconds, no seams or crevices
- No profile key design
- Resistant against most common disinfectants and detergents
- CleanRemind (operating time controlled cleaning/ disinfecting reminder)
- EMC tested according EN 60601-1-2
- Low force and precise feedback key technology, low noise typing
- Compact size, ergonomically low profile design
- Integrated high resolution capacitive touchpad
- CleanLock-function to avoid unintended inputs
- Magnetic fixing option

Switching technology:	Carbon pills on gold contact
Number of Keys:	105 (incl. 12 function keys)
Protection:	IP68 (NEMA 4X)
Key travel:	1,25 mm
Actuation force:	1,2 N
Life:	approx. 2 Mio. actuations
Chassis:	Desktop
Interface:	USB
Power:	20mA@5VDC
Dimensions:	374 x 144 x 15,2 mm
Weight:	0,85 kg
Temperature:	0°C to +60°C
Humidity:	up to 100%
Shock:	3x11 ms pulses (50g) on each of 3 axes
Vibration:	0,04g/Hz, 20Hz-2kHz, 3 hours
Certifications:	CE, FCC Class15, Part B, EN 60601-1-2

Part number	Product name	Features
PS00730	X10-MED-USB-GER	USB, German Layout
PS00731	X10-MED-USB-US	USB, US-Layout

Additional versions and modifications are available on request!

Made by ProKeys e.K.
 Karlsbader Straße 1
 08321 Zschorlau
 Phone +49 (0)3771-5644 356
 Fax +49 (0)3771-5644 357
 www.prokeys.de
 info@prokeys.de