

Silenced GENESI versions without air dryer

MINI BOX 24/5 GENESI S



Code: D-MB075-024-15

Description: Oil free compressor, super silenced, FAD at 5 bar 77 l/min, tank capacity 24 l., without air dryer into mini box kit.

MINI BOX 2016 system warrants the best noise level reduction (49 dB/1m) without any additional connection.

User friendly, compact and light, it can be placed everywhere.

Model								
	Dental chairs	Air dryer	Air Capacity m ³ /h l/min cfm	FAD @ 5bar m ³ /h l/min cfm	Motor HP KW	Volt. V	Freq. Hz	Receiver l GAL
MINI BOX 24/5 GENESI S	1	no	- 130 4,60	- 77 2,7	0,75 0,55	230	50	24 6,34

Model					
	Noise Level dB(A)	Duty cycle	Max pressure bar PSI	Size mm	Net weight Kg lb
MINI BOX 24/5 GENESI S	49	50%	7 101	530x500x800	65 143

