


# Infant Test Lung

## Neonate and infant test lung




Test your neonatal or infant ventilator with the *imtmedical* Infant Test Lung.

### Features

- Two different compliances
- Two different resistances
- Two separate test lungs
- Individual lung selector
- Leak simulator on each lung
- Simple to use

## Specification

- **Tidal Volume:**  
150ml or 300ml
- **Dimensions:**  
120mm x 250mm x 60mm
- **Compliance:** constant over the full range at 0.35 ml/mbar or 0.70 ml/mbar:



## Ordering Information

Partnumber: 500.040.000

TDZ Technika dla zdrowia Sp. z o.o.

ul. Polczyńska 31a  
01-377 Warszawa

tel. mobile:  
+48 500-098-579  
+48 794-098-577

tel.:  
+48 (22) 331-88-29