

MADE IN GERMANY FOR OVER 100 YEARS.
DEDICATED TO DIGNIFIED CARE.

VÖLKER

S 962-2 Vis-a-Vis Hospital Bed





Mobilisation made easy.

Not every patient who is unable or unwilling to leave their bed is completely immobile. Physical complaints or medical interventions often prevent a patient from leaving the bed from the side. But even if this is the case, mobilisation is possible if the bed can be adjusted to allow the patient to sit facing forward so that their feet can make safe contact with the ground straight away. Depending on their capabilities, this can help the patient regain their mobili-

ty and take part in everyday activities. From lying to sitting, through to standing and walking. By successfully combining form and function, Völker hospital beds can provide increased patient comfort while enabling caregivers to work safely and efficiently. The unique Healthcare Design from Völker ensures that all beds fulfil their medical purpose while creating a home-like atmosphere that encourages patient wellbeing and recovery.



→ The bed can easily be adjusted to a sitting position, allowing patients to speak to nurses and visitors at eye level.

Dignified care

Every person has a right to be treated with respect and dignity. Especially people in need of care. We think, develop and manufacture with this in mind.

Patented. The Völker siderail system.

The patented Völker siderail system is characterised by its exceptional versatility. It consists of separate elements that are split on both sides and can be adjusted to either one or two levels. This ensures that patients feel adequately protected from accidentally rolling out of the bed. As the siderails are highly versatile and robust, they can also assist mobilisation.

The side elements are directly integrated into the lying surface. The rails at the head and foot move with the bed as it is adjusted.



Comfortable. The Völker MiS® Micro-stimulation System.

The patented Völker MiS® Micro-stimulation System is based on neuroscience principles. It detects even minor movements made by the patient. The lying surface then feeds these back through the mattress and to the patient's body as tactile stimuli. This continuous feedback increases the patient's body perception, and promotes movement and blood circulation in the skin. Comfort is improved when the patient is sleeping or lying down, mobility is increased and the risk of bed sores is reduced.



Play it safe. With Völker Service.

Völker Service maximises the return on your investment in Völker beds. Our Service team will also take responsibility for a large amount of the work and requirements that arise from the German Medical Devices Act and other legislation.

All of our service technicians are qualified and

experienced testers who are very familiar with the medical devices they inspect. They will test and service your beds meticulously, observing all of the legal requirements.

Any defects are usually fixed immediately and therefore at optimum cost.

Völker Service additionally compiles detailed legal documentation and checks all records to ensure that they are complete.

We would be happy to discuss which Völker Service options would best suit your needs. Simply get in touch.

VÖLKER

MADE IN GERMANY FOR OVER 100 YEARS.
DEDICATED TO DIGNIFIED CARE.

Völker reserves the right to make changes to the design, specifications and models without prior notice. Products may vary slightly from the colours shown here. Subject to printing errors. All information is non-binding.

Not all products/options are available in all countries. For further information on our products or services, please contact your local Völker representative or visit our website.

©2018 Völker | All rights reserved

VÖLKER

www.voelker.de