



P | R | O | D | U | C | T | S

CATALOGUE



BIOLIGHT CO.,LTD.

Add: NO.2 Innovation First Road, Technical Innovation Coast, Hi-tech Zone, Zhuhai, P.R. China Postcode: 519085
Tel: +86-756-3399935 3399999 Fax: +86-756-3399911 E-Mail: Overseas@blt.com.cn

BLT-Y1309-4



Milestones

- 1993 Established as Biolight Electronic Co., Ltd.
- 1996 Introduced ISO9001 Quality Assurance System and certified by TUV.
- 2002 Authorized as “Industrial Basis of National Multi-Parameter Patient Monitor” by National Development & Reform Committee.
- 2004 Patient Monitor products obtained CE certificate.
- 2007 The biggest OEM manufacturer of Patient Monitor in China.
- 2008 Oximeter & Patient Monitor been approved by FDA.
- 2009 New Generation Modular Monitor (Anyview A8) launched.
- 2011 Completes successful IPO & Listed on Shenzhen Stock Exchange.
- 2013 New Q Series Modular Patient Monitor launched.



Patient Monitor





AnyView A8

- 17" TFT LCD ; Resolution: 1280x1024
- Option Module: EMS; MPS; Sidestream / Microstream /Mainstream CO₂ ; Mainstream/Sidestream AG; ICG;C.O.; BIS;CSM;2-IBP; 2-Temp; Nellcor/ Masimo SpO₂
- Global-Leading EMS all-in-one module with display & battery.
- Seamless transport solution, fully meet the clinical demand.
- Touch Screen, DVI & VGA multi-screen display output.



AnyView A6

- 14.1" TFT LCD ; Resolution: 1280x800
- Option Module: EMS; MPS; Sidestream / Microstream /Mainstream CO₂ ; Mainstream/Sidestream AG; ICG;C.O.; BIS;CSM;2-IBP; 2-Temp; Nellcor/ Masimo SpO₂



AnyView A5

- 12.1" TFT LCD ; Resolution: 800x 600
- Option Module: EMS; MPS; Sidestream / Microstream /Mainstream Co₂ ; Mainstream/Sidestream AG;ICG;C.O.; BIS;CSM; 2-IBP; 2-Temp; Nellcor/ Masimo SpO₂



AnyView A3

- 8" TFT LCD ; Resolution: 800x600
- Option Module: MPS; Microstream/Mainstream CO₂ ; Mainstream AG; ICG;C.O.; BIS;CSM;2-IBP; 2-Temp; Nellcor / Masimo SpO₂



AnyView A2E

- 3.5" TFT LCD ; Resolution: 320x240
- Standard Configuration: 3/5 Lead ECG; RESP; SpO₂; PR; 2-Temp; NIBP; Lithium Battery
- Option: 12 lead ECG; 2-IBP; Nellcor / Masimo SpO₂



Emergency Mobile Server (EMS)

- 3.5" TFT LCD display
- Touch Screen & Gravity sensor
- Built-in rechargeable battery
- Seamless transport solution, fully meet the clinical demand.
- Option: MPS (Multi-Parameter Server without display & battery)

NEW

Q7



- 15" TFT LCD, Resolution: 768 x1024
- Module Slot: 3
- Vertical Display for more waveforms
- Fanless cooling system, working quietly.



M9500

- 15" TFT LCD; Resolution: 1024x768
- Standard Configuration: 3/5 Lead ECG; RESP; SpO₂; PR; 2-Temp; NIBP; Lithium Battery
- Option: Sidestream / Microstream /Mainstream CO₂ ; Mainstream/Sidestream AG; 2-IBP; Nellcor/ Masimo SpO₂; Recorder

NEW

Q5/Q3



- 12.1"/10.4" TFT LCD, Resolution: 800x600
- Module Slot: 2
- Plug-and-play design, upgrade easily
- New software, perfect user experience.



M9000A/M8000A

- 12.1" /10.4" TFT LCD; Resolution: 800x600
- Standard Configuration: 3/5 Lead ECG; RESP; SpO₂; PR; 2-Temp; NIBP; Lithium Battery
- Option: Sidestream / Microstream /Mainstream CO₂ ; Mainstream/Sidestream AG; 2-IBP; Nellcor/ Masimo SpO₂; Recorder

NEW

For All Q Series Modular Monitor:




- Standard Configuration: 3/5 Lead ECG; RESP; SpO₂; PR; 2-Temp; NIBP; Lithium Battery; Touch Screen
- Option: 12 Lead ECG; 2-IBP; Sidestream/ Microstream/ Mainstream EtCO₂; Sidestream/ Mainstream AG; Masimo/ Nellcor SpO₂; 2-IBP; 2-Temp; C.O.



M9000/M8000

- 12.1" /10.4" TFT LCD; Resolution: 800x600
- Standard Configuration: 3/5 Lead ECG; RESP; SpO₂; PR; Temp; NIBP; Rechargeable Battery
- Option: Sidestream / Microstream /Mainstream CO₂ ; Mainstream/Sidestream AG; 2-IBP; 2-Temp; Nellcor/ Masimo SpO₂; Recorder



M69

- 12.1" TFT LCD; Resolution: 800x600
- Standard Configuration: 3/5 Lead ECG; RESP; SpO₂; PR; Temp; NIBP; Rechargeable Battery
- Option: Sidestream / Microstream /Mainstream CO₂ ; Mainstream/Sidestream AG; 2-IBP; 2-Temp; Recorder

NEW



M8500

- 8" TFT LCD; Resolution: 800x600
- Standard Configuration: 3/5 Lead ECG; RESP; SpO₂; PR; Temp; NIBP; Lithium Battery
- Option: Microstream/Mainstream CO₂ ; Mainstream/Sidestream AG; 2-Temp; Nellcor SpO₂ / Masimo SpO₂ ;Recorder

M7000

- 7" LCD; Resolution: 234x480
- Standard Configuration: 3/5 Lead ECG; RESP; SpO₂; PR; Temp; NIBP; Rechargeable Battery
- Option: 2-Temp; Recorder



Vital Sign Monitor



V6




- 7" LCD; Resolution: 234x480
- Parameter: BLT Digital SpO₂ ; Nellcor/ Masimo SpO₂ ; PR; Infra-Red ear Temp; NIBP; Microstream/Mainstream CO₂ ; Rechargeable Battery; Recorder



Central Monitoring System

OEM modules

- Multi-Parameter Module Engine
- Digital SpO2 module
- IBP module
- NIBP module
- 3/5/12 lead ECG module
- Thermal recorder





Accessories

- ECG cable and lead wire
- SpO₂ sensor
- NIBP cuff
- Temp probe
- IBP transducers
- EtCO₂ airway adapter & sample lines

M6000C

- Dual screen displays
- Connected to max 66 bedside unit.
- Lan or wireless (WiFi) connected.





E Seires Digital ECG Machine



E80

- 12.1" TFT LCD; Resolution: 800x600
- Size : 400mmx330mmx100mm Weight: ≤7Kg
- Recorder: 12 Channel Folded /Reel thermo printing paper



E70

- 8" TFT LCD; Resolution: 800x600
- Size : 400mmx330mmx146mm Weight: ≤6Kg
- Recorder: 12 Channel Folded /Reel thermo printing paper



E30

- 4.3" TFT LCD; Resolution: 320x480
- Size : 315mmx215mmx77mm Weight: ≤2.25Kg
- Recorder: 3 Channel Reel thermo printing paper



Obstetrical / Gynaecological products



F90/F80 Fetal & Maternal Monitor

- 12.1" TFT LCD; Resolution: 800x600
- Size: 383mmx383mmx124mm
- Parameter: FHR, TOCO, Fetal Movement, Maternal NIBP, SpO₂, PR, ECG, RESP, TEMP.



F50/F30 Fetal & Maternal Monitor

- 7" TFT LCD; Resolution: 800x480
- Size : 328mmx295mmx145mm
- Parameter: FHR, TOCO, Fetal Movement, Maternal NIBP, SpO₂,PR



BLT-870 Digital Colposcope Imaging System

- High resolution Sony imaging system
- Optical Magnification: ≥ 18
- Focus Mode: Auto / Manual Positioning Focus