

# 627/224 Barrier-Free® Examination Chairs

International models



Ritter®  
by MIDMARK



# Midmark Examination Chairs

Accessibility in healthcare goes beyond the ability of patients to access the care they need. The physical accessibility of a facility and the equipment within it plays a significant role. The proper exam chair is important in providing a thorough examination. Midmark Barrier-Free® examination chairs make accessing the exam chair easier, even for those with mobility concerns.



Midmark 627  
Barrier-Free  
Examination Chair

## Standard Features

- 1. Barrier-Free Low Chair Height**  
makes accessing the exam chair easier for all types of patients.
- 2. Powered Back and Height**  
assists in positioning the patient for more comfortable examinations.
- 3. Clean Assist™ Roller System**  
is a retractable roller base that allows you to easily and safely move the chair for cleaning.
- 4. Standard Stirrups**  
are integrated into the exam chair for lower-body examinations.
- 5. Pelvic Tilt**  
makes proper exam positioning easier.
- 6. Front Storage Drawer**  
gives you room to store supplies at the point of care.
- 7. Active Sensing Technology®**  
pauses chair movement when disrupting impact is detected.



**BP Position**  
The Midmark 627 with Patient Support Rails Plus™ help achieve proper patient positioning for blood pressure capture—arm supported at heart height, back supported and feet flat on the floor.



**Elevated Flat Position**  
Easily position patients in a supine or Sims position at a working height that is more comfortable for the provider who prefers to work from a standing position.



**Lithotomy Position**  
Easily position patients for lower-body examinations at a working height that is comfortable for the provider.

# Ritter Examination Chairs



If you're still using a standard, non-adjustable examination table, chances are there have been patients you've had to examine while they sat in a side chair or wheelchair because they couldn't get onto the examination table. This compromises the exam—ensure that your patients receive equal access to quality care regardless of mobility. Ritter examination chairs by Midmark adjust in height to make it easier for even the frailest patient to be seated.

## Help patients help themselves.

A low seat height of 46 cm (18") and wide, flat transfer surface enable patients to transfer onto and off of the exam chair with minimal assistance.



Ritter 224  
Barrier-Free®  
Examination Chair



## Standard Features



### Stirrups and Pelvic Tilt

Help position patients properly and comfortably for lower-body examinations.



## Accessories



### Patient Support Rails Plus™

Helps achieve proper positioning of the patient's arm for blood pressure capture.



### Patient Support Rails™

Provides patients with a 1¼ inch diameter continuous gripping surface for entering, exiting or repositioning on the exam table.



### Assist Arms

Designed to offer patients a stable gripping surface for accessing the chair and a sense of security during positioning.



### Foot Controls

Standard foot control designed to make patient positioning quick, effective and always within reach. (627-011 and 224-011 foot controls shown.)



### Accessory Receivers

Standard receivers are designed to make the integration of accessories easy and efficient—no tools required.



### Articulating Knee Crutches

The unique ball-socket design allows for a maximum range of adjustment, designed to ensure optimum knee and leg support. Available in a non-articulating option.



### Supply Assistant

Provides extra storage for bulky exam supplies. Keeps drapes, gowns and similar items within easy reach.



### Stainless Steel Treatment Pan

Ideal for lower-body procedures or other situations in which removal of fluid and debris is necessary.



### Clean Assist™ Roller System

The retractable roller base allows you to easily and safely move the chair for cleaning.

(Available on the 627-011 model only. Activation levers shown.)



# Specifications

Patient Weight Capacity:  
500 lb (227 kg)

Height:  
Uncompressed - 18" Low/37" High  
(46 cm/94 cm)

Seamless Upholstery Dimensions:  
28" W x 60" L  
(71 cm x 152 cm)

Length with Footrest Extended:  
76" (193 cm)

Pelvic tilt: 5 1/5°

Paper Roll Holder: Holds 24" x 3.56"  
dia. paper roll (61 cm x 9 cm)

Table Weight: 409 lb (185 kg)  
Shipping Weight: 465 lb (211 kg)

## Electrical Ratings:

224-011: 230 VAC, 50/60 Hz, 3.5 A

627-011: 230 VAC, 50/60 Hz, 5 A

## Accessories:

Patient Support Rails Plus™

Patient Support Rails™

Assist Arms

Supply Assistant

Articulating Knee Crutches

Knee Crutches

Stainless Steel Treatment Pan

## 224/627 Comparison Chart

	18" Seat Height (Uncompressed)	500 lb Weight Capacity	Power Base	Standard Rail Accessory Receivers	Active Sensing Technology*	Front Storage Drawer & Pelvic Tilt	Power Back	Clean Assist™ Roller System
224-011	•	•	•	•	•	•		
627-011	•	•	•	•	•	•	•	•

## Midmark & Ritter Standard Upholstery Colors

-  Mist-813
-  Stone-814
-  Restful Path-858
-  Iron Ore-851
-  Lunar Gray-845
-  Obsidian-857
-  Oasis-843
-  Soothing Blue-855
-  Dream-847
-  Mineral-844
-  Healing Waters-848
-  Shaded Garden-853
-  Citrus-849
-  Sandy Retreat-850
-  Curative Copper-846
-  Cranberry-859
-  Dark Linen-856
-  Latte-860
-  Robust Brown-852
-  Deep Earth-854

## Midmark & Ritter UltraFree\* Upholstery Colors

-  UltraFree Stone-865
-  UltraFree Restful Path-868
-  UltraFree Obsidian-867
-  UltraFree Harbor-861
-  UltraFree Wheatgrass-862
-  UltraFree Spice-863
-  UltraFree Cranberry-869
-  UltraFree Dark Linen-866
-  UltraFree Latte-870
-  UltraFree Branch-864

\*ULTRALEATHER® by Ultrafabrics®

### Midmark is an ISO 13485, ISO 9001 and RoHS compliant Certified Company

For more information or a demonstration, contact  
your Midmark dealer or call: 1-800-MIDMARK  
Outside the U.S.A. call: 1-937-526-3662  
or visit our website at midmark.com

Use of ™ and ® in this document denotes Midmark  
trademark rights in the United States. Midmark rights  
in these trademarks may also exist in other jurisdictions.

© 2018 Midmark Corporation  
Midmark Corporation, Dayton, OH.  
Products subject to improvement changes without notice  
Litho in U.S.A. 007-10072-00 Rev. B1 (6/18)

Ritter®  
by MIDMARK



Because we care.