



OXYSTAR

## 4 SEAT

### PORTABLE HYPERBARIC CHAMBER



SIZE:

LENGTH: 215CM

WIDTH: 200CM

WEIGHT:

563KG

POWER:

1200W

WORKING AIR PRESSURE:

7PSI (1.5ATA)

COMPRESSOR AIR FLOW:

133L/MIN

O2 LPM:

10L/MIN

TWO COMPRESSOR

ONE CONCENTRATOR

FOUR OUTPUT VALVES

OXYGEN TUBE:


92%~98%

DIFFUSE OXYGEN SUPPLY:

28%~35%

NEGATIVE IONS NUMBER:

100000PCS/CM3



HBOT POLSKA sp. z o. o.

+48 518 412 200

infohbot@gmail.com