




OXYSTAR

## 8 SEAT

### PORTABLE HYPERBARIC CHAMBER



SIZE:

LENGTH: 255CM

WIDTH: 213CM

WEIGHT:

739KG

POWER:

1650W

WORKING AIR PRESSURE:

7PSI (1.5ATA)

COMPRESSOR AIR FLOW:

133L/MIN

O2 LPM:

10L/MIN

TWO COMPRESSOR

TWO CONCENTROR

FOUR OUTPUT VALVES

OXYGEN TUBE:

92%~98%

DIFFUSE OXYGEN SUPPLY:

28%~35%

NEGATIVE IONS NUMBER:

100000PCS/CM3

**HBOT POLSKA sp. z o. o.**

**+48 518 412 200**

**infohbot@gmail.com**