

FlowLab

Software for CITREX and FlowAnalyser

analyser
the art of measuring

Reporting and documentation software

View your flow, pressure and volume measurements comfortably on your computer monitor while simultaneously saving respiratory parameters. FlowLab displays your measurement results graphically in real-time curves or numerically. Trending reports of up to 100 hours can be generated to verify long-term ventilator functionality. To facilitate your device management, test reports can also be created and saved electronically or printed.

« Software package for real-time curves,
100h trending and customized reporting. »



Test Report

- Test report printouts with one click
- User defined configuration
- Logo insertion available
- FlowAnalyser data automatically retrieved
- Various input options for each tested ventilator
- Unique control number for each report

Panels

- Displayed in time-relation or as a loop
- Various cursors measure the curves
- Unique trigger used to display real-time curves in Single Shot, Norm or Auto mode
- User defined layout and colours
- Option of setting a title, printing and saving
- Simultaneous display of up to 6 different curves

Trending

- Up to 100 hours of data logging
- User defined trending interval
- Up to 10 values simultaneously
- Data export to Excel, etc.
- User defined layout and colours
- Option of setting a title, printing and saving
- Automatic file size calculation

Numerics

- All measurements displayed on one page or combined with panels
- User defined layout and colours
- Statistical data including mean, max and min for each value
- Input of target value with tolerances
- Automatic verification of Pass/Fail
- Up to 20 values can be displayed simultaneously

Languages

- Available in 15 languages
- All languages included – free of charge
- Languages from every continent

