

**NEW
INNOVATION**

Geratherm®

UniqueTemp° **X3**

X3 Applications – one system

Warming blankets and mattress




- Up to 4 OR warming blankets plus 1 OR table mattress
- Reusable

Control unit



- Intuitive menu navigation via touch-screen

Fluid warming



- Warming tube for standard infusion sets
- Reusable



Temperature sensor

- Ear sensor for precise measuring of core temperature
- Non-invasive use

Geratherm[®]

UniqueTemp[°] 3x

Technical facts	Geratherm UniqueTemp [°] 3x						
control unit	digital control for all possible combinations						
Art. no.	8001340800						
Operation	Touch-screen menu navigation and selection						
Firmware	updateable to the latest version						
Dimensions weight approx.	H 23.0 cm x B 11.0 cm x T 16.5 cm 3.6 kg						
Voltage	Input voltage: 100 - 240 V, AC Output voltage: 35.5V, DC						
Safety features	Autocheck, temperature measurement between warming system and patient, constant temperature monitoring, autostop in cases of overheating						
Connections to the blankets	Push & Pull quick-release coupling (as bus-system)						
Warming segments	Arm/Shoulder blanket	Torso blanket	Leg blanket		OR table mattress	recovery blanket	
			right	left		right	left
Design							
Art. no. blanket	8001340202	8001340301	8001340408	8001340409	8001340600	8001350500	8001550501
Art. no. cover	8001340203	8001340303	8001340410	8001340411	8001340601	8001350502	8001350503
Temperature pre-setting	30 °C - 40 °C; in steps of 0.1 K; separate settings for each blanket				30 °C - 39 °C in steps of 0.1 K	30 °C - 40 °C; in steps of 0.1 K; separate settings for each blanket	
Temperature display	Desired and actual temperature visible on colour display (TFT 3.5")						
Material	Special textile, blood and fluid resistant						
Cleaning	wipeable, resistant to disinfectant, covers in addition removeable and washable up to 95 °C						
Weight approx.	blanket: 0.9 kg cover: 0.4 kg	blanket: 0.75 kg cover: 0.2 kg	blanket: 0.75 kg cover: 0.2 kg		mattress: 2.2kg cover: 0.2kg	blanket: 2.3kg cover: 0.9kg	
Dimensions approx.	175 cm x 32 cm	45 cm x 92 cm	45 cm x 92 cm		93 x 49 cm	120 x 200 cm	
warming tube	actively heated profile for fluid warming to be used with standard infusion sets (5 mm)						
Art. no.	8001340700						
temperature settings	33 °C – 40 °C						
Temperature display	Desired and actual temperature visible on colour display (TFT 3.5")						
Material	silicone						
Cleaning	wipeable, resistant to disinfectant						
Weight approx.	0.35kg						
Dimensions approx.	150cm						
temperature sensor	sensor for measuring the patient's temperature in the outer ear canal						
Art. no.	8001340700						
measuring range	0 °C - 45 °C (accuracy: 35-42 °C: 0.15 K; 0-34.9 and 42.1 - 45 °C: 0.2 K)						
Material	Silicone, elastomer						
Cleaning	wipeable, resistant to disinfectant; hygienic disposable caps						

Subject to technical modifications