


DynaTrac™ Retraction System

FOR OPTIMAL SURGICAL EXPOSURE AND STABILIZATION



- Precise, dynamic retraction for superficial and deeper retraction
- Retractor frame offers versatility for a variety of surgical procedures
- Cost-effective, reusable retractor frame and blades
- Simplifies dissection and overall surgery
- DynaTrac™ frames are lightweight, blue anodized aluminum
- FlexArm™ Plus Holding & Positioning System affords maximum stabilization, particularly during procedures in the pelvic region

EFFICIENCY • STABILITY • VERSATILITY

DynaTrac™ Retraction Systems

The Mediflex® DynaTrac™ Retraction System enables precise retraction for optimal surgical exposure and access. The versatile, lightweight DynaTrac™ Retractor Frame permits use of elastic Surgical Stays for superficial retraction as well as retractor blades with ratchets for deeper retraction.

DynaTrac™ Frames in 3 configurations to accommodate multi-specialty procedures

Horseshoe


Ref. No. 78000-HP



For Gynecology and Urology

Square

Ref. No. 78001-HP



For Colo-Rectal, Urology, Gynecology, General Surgery and Head & Neck

Figure 8

Ref. No. 78002-HP



For Gynecology, Urology and Head & Neck

Disposable Elastic Surgical Stays



Ref. No. 4005S



Ref. No. 4005B



Ref. No. 4012B

The FlexArm™ Plus Holding & Positioning System provides maximum stability of the DynaTrac™ Retraction System, and can additionally be used for scope or instrument holding during laparoscopic procedures



Ref. No. 99054-QCLR



Shown with Ref. No. 72162

ORDERING INFORMATION

PRODUCT	REF. NO.	DESCRIPTION
DYNATRAC™ RETRACTOR FRAMES with HEX POST	78000-HP	Horseshoe Frame, 8.4" (21.3cm) x 7.8" (19.8cm)
	78001-HP	Square Frame, 7.4" (18.8cm) x 7.4" (18.8cm)
	78002-HP	Figure 8 Frame, 15.3" (38.9cm) x 8.9" (22.6cm)
DISPOSABLE ELASTIC SURGICAL STAYS	4005S	5mm Sharp Hook, Sterile, 4 Stays/Pk, 5 Pks/Bx
	4005B	5mm Blunt Hook, Sterile, 4 Stays/Pk, 5 Pks/Bx
	4012B	12mm Blunt Hook, Sterile, 4 Stays/Pk, 5 Pks/Bx
BLADE RATCHET	72294	Mini-Bookler® Tilt Ratchet
	72256	Double Hook, Blunt, .2 x .4" (.5 x 1cm)
MINI-BOOKLER® HOOK RETRACTORS	72252	Triple Hook, Blunt, .35 x .4" (.9 x 1cm)
	72257	Quad Hook, Blunt, .5 x .35" (1.3 x .9cm)
	72150	Quad Hook, Deep, Blunt, .7 x .55" (1.8 x 1.4cm)
	72150-S	Quad Hook, Deep, Sharp, .7 x .55" (1.8 x 1.4cm)
	72276	Deaver, ½ x 2" (1.3 x 5.1cm)
	75638	Deaver, 1 x 3½" (2.5 x 8.9cm)
	75639	Deaver, 1 x 5" (2.5 x 12.7cm)
MINI-BOOKLER® SURGICAL RETRACTOR BLADES	72263	Kelly, ½ x 2" (1.3 x 5.1cm)
	72269	Kelly, 1 x 1½" (2.5 x 3.8cm)
	72271	Kelly, 1½ x 2" (3.8 x 5.1cm)
	72259	Malleable, ½ x 2½" (1.3 x 6.4cm)
	72273	Malleable, ½ x 4" (1.3 x 10.2cm)
	72272	Malleable, ¾ x 4" (1.9 x 10.2cm)
	72274	Malleable, 1 x 4" (2.5 x 10.2cm)
	75642	Malleable, 1 x 5" (2.5 x 12.7cm)
	75633	Rectal, ⅞ x 2⅞" (2.2 x 7.3cm)
	75634	Rectal, ⅞ x 3⅞" (2.2 x 9.8cm)
	75635	Renal, Wiley, 1 x 7" (2.5 x 17.8cm)
	75636	Vaginal, Lateral, 1 x 3" (2.5 x 7.6cm)
	75637	Vaginal, Posterior, 1 x 4½" (2.5 x 11.5cm)
FLEXARM™ PLUS SYSTEM	99054-QCLR	FlexArm™ Long Reach System, QC Hex Fitting
HOLDING ATTACHMENTS	72162	10mm Quick-Grip Automatic Holding Tip
	72163	5mm Quick-Grip Automatic Holding Tip

Mediflex®
SURGICAL PRODUCTS

250 Gibbs Road • Islandia, NY 11749

© 2019 Mediflex® Surgical Products, Islandia, NY. All rights reserved.

Mediflex® is a registered trademark of Mediflex® Surgical Products, a division of Flexbar Machine Corp.

www.mediflex.com

USA 800.764.0474

+1.631.582.6424

info@mediflex.com

DRS-0919