




2020

Manipulatory macity



Uterusmanipulator, adjustable 90°

Uterusmanipulator, rigid

For a gentle surgery

Colpotomy set

Uncomplicated use in gynecology



Atraumatic tips






Manipulatory macicy „sztywne” bez regulacji kąta uniesienia końcówki:

The new Uterusmanipulator, rigid

For a gentle surgery

The rigid manipulator serves a simplified manipulation of the uterus at the hysterectomy. The different caps (Ceramic, Teflon and PEEK) can all be used in the high frequency surgery.



Accessories: 02-00052 Closing sponge, sterile, disposable, pack of 5
02-00053 Adapter, 11 parts, WL 20 mm - 120 mm

02-00045-12 Ceramic cap, Ø 32 mm
02-00045-13 Ceramic cap, Ø 36 mm
02-00045-14 Ceramic cap, Ø 42 mm



NEW
for use with
ultrasound

02-00045-24 Teflon cap, Ø 32 mm
02-00045-25 Teflon cap, Ø 36 mm
02-00045-26 Teflon cap, Ø 42 mm



02-00045-23 PEEK cap, large






Manipulatory macicy z regulacją kąta uniesienia końcówki:

Hysterectomy

The new Uterusmanipulator, adjustable 90°

For a gentle surgery

The new Uterusmanipulator for laparoscopic total or subtotal Hysterectomy.



Accessories: **02-00052** Closing sponge, sterile, disposable, pack of 5
02-00053 Adapter, 11 parts, WL 20 mm - 120 mm

02-00045-12 Ceramic cap, Ø 32 mm
02-00045-13 Ceramic cap, Ø 36 mm
02-00045-14 Ceramic cap, Ø 42 mm



NEW
for use with
ultrasound

02-00045-24 Teflon cap, Ø 32 mm
02-00045-25 Teflon cap, Ø 36 mm
02-00045-26 Teflon cap, Ø 42 mm



02-00045-23 PEEK cap, large





Good visualisation of vagina and rectum.


The vaginal/rectal probe serves the pelvic floor reconstruction and has to be inserted vaginally or rectally. For the attachment of mesh implants it allows a better visualisation of vagina and rectum.



Vaginal vault prolapse before operation

Treatment with vaginal/rectal probe and mesh

Atraumatic tips



— Ordering information —

02-00093 Vaginal probe Ø 40 mm x 35,5 mm, reusable

02-00092 Rectal probe Ø 28 mm x 34,5 mm, reusable

Stanisław Wawrzyniak

Tel. +48 883 013 093

s.wawrzyniak@yomed.pl

Fax: +48 32 352 09 33

YOMED Joanna Wawrzyniak

ul. Poleska 18/2, 40-733 Katowice,

NIP 5860055782, biuro@yomed.pl, www.yomed.pl