


Novel Coronavirus (2019-nCoV) IgG/IgM Test Kit

Colloidal Gold





COVID-19 can be screened by measuring the antibody IgG & IgM for novel coronavirus. Antibody IgM appears in the early stage of infection, while Antibody IgG in the late stage.

Easy

Only three steps to read results

Efficient



Reliable results in less than 15 minutes



Novel Coronavirus (2019-nCoV) IgG/IgM Test Kit (Colloidal gold)

A Diagnostic Solution for COVID-19

Interpretation



Packaging

IgG and IgM results in one test
25 tests per kit
Each test card in a pouch
2.5mL diluent
Storage at room temperature

