

**DIMENSIONS (±5%)**

Mattress	1970 x 860 mm
Overall dimensions	2180 x 1000 mm
Minimum/Maximum height	240 / 780 mm
Platform/Side Rails height	373 mm

\* Casters Ø100mm on ABS mattress platform  
HPL mattress platform: 220 / 760 mm

**ARTICULATIONS (±5%)**


Back-rest maximum angle	70°
Cardiac Chair	88°
Leg-rest maximum angle	32°
Foot-rest maximum angle	28°
Maximum tilt	
Trendelenburg	18°
Anti-Trendelenburg	18°

**CAPACITIES (±5%)**

Patient maximum load	185 Kg
Safe Working Load (SWL)	250 Kg
Overall bed weight	97 Kg
Electric tension	100/240VAC-50/60Hz
Maximum power	200 W
Batteries maximum power	18 W

**REGISTRY AND CONFORMITY**

Officially registered as Medical Devices Class I producer over the No. 125/DM/2017/01/217.  
CE Mark.  
Product follows UE directive 2007/47/EE, annex IX.  
Quality Management certification EN ISO 9001:2008  
Medical Devices certification EN ISO 13485:2016



indicative colors only, check with real samples.

IMO-INDUSTRIAS METALURGICAS, S.A.  
www.imo.com.pt

IMO PORTO  
Rua Adolfo Casais Monteiro, 117  
4050-520 Porto  
Tel.: 222 009 767  
Fax: 222 000 647  
porto@imo.com.pt

HEADQUARTERS & FACTORY  
Rua Dr. Aurélio Teixeira de Sousa, 575  
4650-312 Rande  
PORTUGAL  
T. 255 340 220  
F. 255 340 239  
sede@imo.com.pt

IMO LISBOA  
Av. Defensores de Chaves, 87-A  
1000-115 Lisboa  
Tel.: 217 963 478  
Fax.: 217 975 822  
lisboa@imo.com.pt

AUTHORISED DISTRIBUTOR

COPYRIGHTS - IMO, S.A. | CATALOGUE 02.018.3 - JUNE 2018 | WE RESERVE THE RIGHT TO CHANGE PRODUCTS SHOWN WITHOUT PREVIOUS NOTICE



medical beds worldwide


**MATRIX U24**  
highlights




Trendelenburg / Anti-Tredelenburg position



Cardiac Chair position



8 functions portable control.  
Backlit.




Supervisor control

# MATRIX U24

**SAFETY SIDE-RAILS**  
collapsible metal side-rails  
(EN 60601-2-52 standard)

**MATTRESS PLATFORM**  
removable ABS modules



**WHEELS CLUSTER**  
4 casters Ø100mm with central lock

ultra  
24 cm  
low