OxySmart

Bluetooth Fingertip Pulse Oximeter









- Wireless monitoring & data storage on APP
- Artifact removal and anti-motion
- Accurate and sensitive

Features

- OLED display of SpO2, PR, PI, Pulse bar and waveform
- Special design with splash-proof and drop-resistant
- Spot check and continuous measuring mode
- Audible and visual alarm with low battery indication

Specifications

SpO2		Power supply
Transducer	Dual-wavelength LED sensor	2 x AAA alkaline battery
Range	35%~100%	
Accuracy	≤ 2% from 70% to 100%	Perfusion Index (PI)
		Display range 0%~20%
Pulse Rate		Data Storage
Range	30bpm~240bpm	Up to 12 groups
Accuracy	±2bpm or ±2%	<u> </u>
		Measuring Mode
Default Alarm Limit		30s spot check and continuous mode
SpO2 low limit	90%	
Pulse rate	High: 120bpm	Pulse Rate Analysis(Spot check mode)
	Low: 50bpm	12 types preset analysis



Shenzhen Viatom Technology Co., Ltd.

www.viatomtech.com

marketing@viatomtech.com