



# gMotio



# gMotio

The gMotio gynaecology couch was designed following 2 years of study with health professionals in Europe and the Middle East. It has been tested by around 20 doctors and 100 patients. gMotio has been designed to provide effective and comfortable treatment, both for the practitioner and for their patients. With its extensive experience, Promotal has created a synergy between two featured products: the Elite gynaecology couch and the eMotio examination couch. gMotio is already considered a must-have gynaecology couch.





## EQUALITY OF TREATMENT

### MAKE LIFE EASIER FOR YOUR PATIENTS WITH REDUCED MOBILITY WITH:

- > Starts at 57cm from the ground
- > Electrical support of all patient movements Electrical
- > moving footrests (depending on model)



## EFFICIENCY

### IMPROVE THE FLOW OF YOUR EXAMINATION WITH:

- > An integrated leg rest extension for a trendelenburg position (opt)
- > 3 programmable positions (9 as accessories)
- > Autonomous movements: Press once to change positions
- > Intuitive control and its already configured patient start position
- > Integrated lighting in the table (depending on model)
- > Standard seat rails



## ERGONOMICS

### ALWAYS PROVIDING THE BEST POSITION FOR YOU AND YOUR PATIENT:

- > Precise work on the seat design and foam density
- > Variable height from 57 to 118,5cm (with tilt)
- > Precise and adjustable tilt, easy speculum placement
- > Controls at the right place on the backrest and on the feet
- > Built in retractable wheels for short trips



## SAFETY

### EASY CLEANING OF YOUR COUCH:

- > Seamless upholstery
- > Removable seat and backrest for easy cleaning
- > Cataphoresis treatment of the structure for longevity and hygiene of the product
- > Removable ABS tray
- > Electrically activated wheels by sustained contact

### SAFEGUARD YOUR EXAMINATIONS:

- > Safety sensors
- > Safety battery
- > Dual control

# gMotio



**ELECTRICAL FOOT RESTS  
+ INTEGRATED LIGHTING**  
— Ref 15640-01 (pictured above)

**ELECTRICAL FOOT RESTS**  
— Ref 15630-01



**STIRRUPS ON CLAMPS**  
— Ref 15620-01



**LEG RESTS ON CLAMPS**  
— Ref 15626-01



**GOEPEL ON CLAMPS**  
— Ref 15624-01

# gMotio

- OPTIONS
- ACCESSORIES



Integrated mains socket and equipotential terminal  
(EUROPE, UK, US) Ref. 40400-1  
(DK, IT, AUS, CH) Ref. 40400-2



Urology tray  
(replacing the std tray)  
Ref. 2061-30



Stainless steel tray  
(replacing the std tray)  
Ref. 15966-01



Leg rest extension  
Ref. 15672-01



Manual remote control  
(9 programmable positions)  
Ref. 15031-10



Pair of coated stirrups  
(add x2 clamps 878-12)  
Ref. 12850 (+ 2 x 878-12)



Pair of plastic leg rests  
(add x2 clamps 878-12)  
Ref. 840 (+ 2 x 878-12)



Pair of covers for leg rests  
Ref. 861



Pair of goepfels  
(add x2 clamps 878-12)  
Ref. 845  
(+ 2 x 878-12)



Single armrest  
(add x1 clamp 878-12)  
Ref. 2536  
(+ 878-12)



Armrest cover  
(for model 2536)  
Ref. 2537



Footrest hoop cover  
(for models 15630-01 & 15640-01)  
Ref. 15861-01



Adjustable headrest  
(iQuest 3055-20 and gmotio)  
Ref. 30679-01



PVC footrest protection  
Ref. 15905-01



PVC seat protection  
Ref. 15906-01



PVC extension leg rest protection  
Ref. 15907-01



Self-locking serum holder  
rod 2 hooks  
Ref. 2985-01 (+ 878-12)



Serum rod holder  
2 hooks with stop  
Ref. 985-01 (+ 878-12)



Splash guard  
Ref. 12677-01



Colpo support. Zeiss d.38  
for gMotio\*  
Ref. 15840-01

\*suppresses autonomous movements



Colpo support. Pivoting  
Leisegang for gMotio\*  
Ref. 15810-01

\*suppresses autonomous movements



Olympus right colpo  
support\* (elite/elansa 11614-01)  
Ref. 15850-01

\*suppresses autonomous movements



Olympus left colpo  
support\* (elite/elansa 11614-01)  
Ref. 15851-01

\*suppresses autonomous movements



Stainless steel backrest  
rail (pair) 25x10 mm  
Ref. 15879-01



Sliding clamp for rail 25x10  
mm, ø 16/18 mm  
Ref. 878-12



Rotary clamp for rail  
25x10 mm  
Ref. 879-10






Lamp holder for clamps  
Ref. 957 (+ 878-12)




Lamp  
(Refer to the documentation)

# gMotio



15620-01 / 15624-01 / 15626-01



15630-01 / 15640-01

## UPHOLSTERY COLOURS

32 - BLUE		CLASSIC
93 - COTTAGE		
35 - RED		
34 - YELLOW		
91 - WHITE		
92 - MARINE		VITAMIN
56 - DARKNESS		
57 - ANTISTATIC		
64 - GRAPE		
67 - LAGOON		
66 - SWEET		COSY
63 - CEDRAT		
80 - SKY		
81 - SEA GREEN		
65 - LILAC		
96 - TIN		
82 - MOUSE GREY		
83 - NATURAL		
84 - CHOCOLATE		