



YOU MEASURE

NeuroLight[®]

The new era pupillometer for neurological monitoring



SIMPLE AND PRECISE ROUTINE CHECKS
FOR PUPIL SIZE AND REACTIVITY

- Quantitative measurements
- Non invasive
- Time saving

idmed
an eye on your patient

Unique features

+ Ergonomic design
Easy to operate

+ Large touch screen
Intuitive use

+ Long lasting battery
1 week of autonomy

+ Integrated barcode reader
Simplified patient identification

+ Data transfer
Save time

+ Opaque eyecup to ambient light
Reliable and reproducible measurements

+ High definition camera
Precision 0,1mm



FL Calibrated light stimulation

Graph Pupil reactivity

0-24h 24h trend chart

L/R Anisocoria detection

Person + 10 000 measurements

NeuroLight® and IDMED® are registered trademarks owned IDMED SAS.
Photos: ©AdobeStock - ©IDMED - Copyright 2021 - IDMED SAS. All rights reserved.

Ref: DC-NL-V5-EN - Revision: 03/2021

idmed
@IdmedNews · in

45 rue Joliot Curie - Hôtel Technologique
13013 Marseille FRANCE
Tel. : +33(0)4 91 11 87 84 • Email : info@idmed.fr