

A lightweight, wieldy, clean free solution for Intubation

- Lightweight
- Smallest Diameter
- Disposable Sheath with Patent Design
- Rotatable Display
- Versatile Video Output



Slim Stylet Video Scope

StyletViu® provides adult size, child size and infant size video stylet, the first of it's kind in the world. It generates real-time video images and most importantly, it is safe, reliable, and easy to use for endotracheal intubation in any difficult situations.

Features

- Lightweight with world's lightest 80g video scope, easy to hold.
- Thin 2.4/4.5mm for ET tube 2.6, 3.5, 7.0
- Safe design with camera inside stylet.
- Micro camera has 2/2 LEDs in the distal end with wide 90° coverage.
- Rotatable 2.4" vertical view display with angle adjust by thumb.

Technical Specification



Model Number	CA76L-TF24-L220	/ CA76L-TF45-L350
Disposable Sheath	3.0 mm with patent design	/ 6.0 mm with patent design
Sensor	Miniature EndoCam™ SD sensor	
Field of View	90°	
Distal End	2.4	/ 4.5 mm
Working Length	220	/ 350 mm
Type	Reusable	/ Reusable
Display System	2.4" Vertical IPS screen plug-on control system	
Resolution	SD	