

## FlexArm<sup>™</sup> Plus System for Whipple / Bariatric Procedures

Mediflex offers a holding & positioning system to provide rigid, low-profile liver retraction for optimal surgical field access - without robotic interference



FlexArm™ Plus goes where human arms can't...



- FlexArm<sup>™</sup> Plus permits low-profile positioning without blocking or hindering robot arm movement
- Lapro-Flex<sup>®</sup> articulating retractor 'pushes' the left liver lobe to optimize surgical field view
- Eliminates robotic arm or staff requirement to hold retraction
- Cost-effective, reduces OR time

ORDERING INFORMATION		
Product	Ref. No.	Description
FLEXARM SYSTEM	99054-QCLR	FlexArm Plus with Quick-Connect Hex Fitting - Long Reach Post and OR Rail Clamp
HOLDING TIP	69706	Stainless Steel Adjustable Tip Holding Capacity: 3.5mm – 15mm
LAPRO-FLEX <sup>®</sup> RETRACTORS	91682-A	60mm Angled Triangle Base, 5mm x 32cm Shaft Length
	91683-A	80mm Angled Triangle Base, 5mm x 32cm Shaft Length

Lapro-Flex<sup>®</sup> side profile of angled base

## Mediflex Surgical Products

Toll-free (USA): 1.800.764.0474 Phone: +1.631.582.6424 Fax: +1.631.582.4520 www.mediflex.com

## NO INTERFERENCE • VERSATILE • FLEXIBLE

## FAS/WB-0518